

## Product datasheet for **MG201358**

### Ccdc12 (NM\_028312) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ccdc12 (NM\_028312) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Ccdc12  
**Synonyms:** 2700094L05Rik; C76605  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG201358 representing NM\_028312  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCGGCTGCTCCTGCCGGTGTGGCCGCTAGAGGAAGAGGCGCTGCGCGGAAAGAACGGCTGAAGG  
 CCCTCCGGGAGAAAACCGGGCGCAAGGACAGGGAAGATGGGGAGCCACAGACCAAGCAGCTCAGAGAAGA  
 GGGAGAGGAAGTTGGGAAGCACAGGGGACTCAGGCTCCGGAATATGTCCCAGAGGATGAAGACCTGAAG  
 AGGAGGGGGTCCCCAGGCCAAACCAGTGGCAGTGGAGGAGAAGGTGAAGGAGCAGCTGGAGGCAGCCA  
 AGCCAGAGCCTGTCATTGAGGAAGTGGACCTGGCCAACCTGGCACCCGAAACCTGATTGGGATCTCAA  
 GAGAGATGTGGCCAAGAAGCTAGAGAAGCTAGAAAAGCGGACCCAGAGAGCCATCGCGGAGCTGATCCGT  
 GAGCGGCTAAAAGGCCAGGAGGACAGCCTGCCCTCTGCAGTGGATGCCACCACCGGACAAGAGGCTGTG  
 ACTCTGAC

**ACGCGT**ACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG201358 representing NM\_028312  
 Red=Cloning site Green=Tags(s)

MAAAPAGVGRLEEEALRRKERLKERLREKTGRKDREDGEPQTKQLREEGEEVKGHRGLRLRNYVPEDEDLK  
 RRRVPAKPVAVEEKVKEQLEAAKPEPVIEEVDLANLAPRKPDPDLKRDVAKKLEKLEKRTQRAIAELIR  
 ERLKGQEDSLPSAVDATTGQEACDSD

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_028312

**ORF Size:** 498 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_028312.1](#), [NP\\_082588.1](#)

**RefSeq Size:** 783 bp


**RefSeq ORF:** 501 bp

**Locus ID:** 72654

**UniProt ID:** [Q8R344](#)

Cytogenetics: 9 F2

Product images:



Circular map for MG201358