

Product datasheet for **MG201423**

Dctpp1 (NM_023203) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dctpp1 (NM_023203) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Dctpp1
Synonyms: 2410015N17Rik; AI854235; RS21-C6
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201423 representing NM_023203
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCCACGGCTGGTGACGGTGAGCGGTAAGTGTGGGACAAGAAGACTCTGCTGCTGCCGCCGTTCA
 GATTCAGTCCGGAGCCCACGCTTGGAGACATCCGACGCTCCATGCAGAGTTGCAGCTGAGCGGGACTG
 GGAGCAGTTCCACCAGCCTCGGAACCTGCTGCTAGCCTTGGTGGGAGAAGTGGCGGAGCTGGCAGAACTC
 TTCCAGTGAAGTCTGATACAGAGCCTGGTCTCAAGCATGGCCACCAAAGGAGAGAGCAGCCCTTCAGG
 AGGAGCTTAGCGACGTCCTCATCTACCTGGTAGCATTAGCAGCCGCTGCCATGTGGATTACCTCAAGC
 AGTAATCTCCAAAATGGACACCAACCGCAACGTTACCCAGTCCACCTGTCCCGAGGTTCTGCCTGCAAG
 TACTACTGACTTGCCCGTGGGACCATCTCTGAAAACAGGCTGTGGGGGCTGGAGACCCTGCCTCGGAGC
 TGAGAGACCAGGCTTCCACA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201423 representing NM_023203
 Red=Cloning site Green=Tags(s)

MSTAGDGERGTVGQEDSAAARPFRRFSPEPTLEDIRRLHAEFAAERDWEQFHQPRNLLLALVGEVGEAEL
 FQWKSDETPGPQAWPPKERAALQEELSDVLIYLVALAARCHVDLPQAVISKMDTNRQRYPVHLSRGSACK
 YTDLPRTISENQAVGAGDPASELRDQAST

TRTRPLE - GFP Tag - V


Restriction Sites: Sgfl-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM_023203

ORF Size: 510 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_023203.1, NP_075692.1</u>
RefSeq Size:	648 bp
RefSeq ORF:	513 bp
Locus ID:	66422
UniProt ID:	<u>Q9QY93</u>
Cytogenetics:	7 F3
Gene Summary:	Hydrolyzes deoxynucleoside triphosphates (dNTPs) to the corresponding nucleoside monophosphates. Has a strong preference for dCTP and its analogs including 5-iodo-dCTP and 5-methyl-dCTP for which it may even have a higher efficiency. May protect DNA or RNA against the incorporation of these genotoxic nucleotide analogs through their catabolism. [UniProtKB/Swiss-Prot Function]