

## Product datasheet for **MG201436**

### Zfand2a (NM\_133349) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Zfand2a (NM\_133349) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Zfand2a  
**Synonyms:** AA407930; Airap; AU016206  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG201436 representing NM\_133349  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGAGTTTCCTGACTTGGGGAAGCACTGTTTCAGAACCACTTGCAAACAGCTAGATTTTCTGCCAATAA  
 CATGCGACGCCTGTAACAAGATTTCTGTAAGACCATTTTTCTACGTGGGTATAAGTGTCCCTTTGC  
 GTTCAAGAAGGATGTGCAGGTGCCTGTGTCTCTCTGCAATGCCCCATCCCTGTGAAGAGAGGTGAG  
 ATCCCAGATGTGGTGGTGGCGAGCACATGGACAGAGATTGCACCTTTCACCCTGGGAGGAACAGAAACA  
 AGTTTTTCACACATCGCTGCTCCAAAGAGGGATGCAGGAAGAAGGAGATGCTGCAGCTGGCGTGTGCTCA  
 GTGCCATGGCAACTTCTGCATTCAGCACAGGCATCCTCTGGACCAATTGCCAGGCTGGGAGCAGCTCA  
 GCCAGCAGAGGGAGGACTTCTACATCCAGAGCTGCTGAGCAGAAGCCATCGGGAGTCAGTTGGCTGGCC  
 AGCGGCTCAGGCGCACAGTGAAG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG201436 representing NM\_133349  
 Red=Cloning site Green=Tags(s)

MEFPDLGKHCSEPTCKQLDFLPITCDACKQDFCKDHFVYVGHKCPFAFKKDVQVPVCPNAPIPVKRG  
 IPDVVVGEHMDRDCTFHPGRNRNKVFTHRCSKEGCRKKEMLQLACAQCHGNFCIQHRHPLDHNCQAGSS  
 ASRGRSTSRAAEQKPSGVSFLAQLRRTVK

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)



<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_133349.3</a></u> , <u><a href="#">NP_579927.1</a></u>
<b>RefSeq Size:</b>	3115 bp
<b>RefSeq ORF:</b>	516 bp
<b>Locus ID:</b>	100494
<b>UniProt ID:</b>	<u><a href="#">Q9JII7</a></u>
<b>Cytogenetics:</b>	5 G2