

Product datasheet for **MG201470**

Rps16 (BC090618) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rps16 (BC090618) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Rps16
Synonyms: AA420385; AI317031
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201470 representing BC090618
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTGGAGCGGCCGCTGAGAGACCTTTCCGGTCGCGGCGCTGCGCGGTGTGGAGCTCGTGCTTGTGC
TCGGAGCTATGCCGTCCAAGGTCGCTGCAGTCCGTGCAGGTCTTCGGACGCAAGAAAACAGCCACGGC
TGTGGCCACTGCAAACGGGAAATGGGCTCATCAAGGTGAACGGACGTCCCCTGGAGATGATCGAGCCG
CGCACGCTGCAGTACAAGTTACTGGAGCCTGTTTCTGCTTCTGGCAAGGAGCGATTTGCTGGTGTGGATA
TTCGGTCCGTGTGAAGGTGGTGGACATGTGGCCAAATTTATGCCATCCGACAGTCCATCTCAAAGGC
CCTGGTAGCTTATTACAAAAATATGTGGATGAAGCCTCCAAGAAGGAGATCAAAGATATCCTCATCCAA
TACGATCGGACCCTGCTTAGCTGACCCCGTCGCTGCGAATCCAAAAGTTTGGAGGTCTGGTGCC
GTGCCGATACCAGAAATCCTACCGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201470 representing BC090618
Red=Cloning site Green=Tags(s)

MLERPLREPFVVAALRGVELVVLVGLAMPKSGPLQSVQVFRKKTATAVAHCKRGNLIKVNRPLEIIEP
RTLQYKLLLEPVLLLKGERFAGVDIRVRVKGGHVAQIYAIRQSISKALVAYYQKYVDEASKKEIKDILIQ
YDRTLLVADPRRCESKKFGGPGARARYQKSYR

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



Cloning Scheme:


ACCN: BC090618

ORF Size: 518 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC090618](#), [AAH90618](#)


RefSeq Size: 600 bp

RefSeq ORF: 518 bp

Locus ID: 20055

Cytogenetics: 7 A3

Product images:



Circular map for MG201470