

Product datasheet for **MG201509**

Wrb (NM_207301) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Wrb (NM_207301) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Wrb
Synonyms: 5530402J05Rik; C030018G21Rik; Chd5
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201509 representing NM_207301
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGCGCGTCCGAGACGGACCGCTGGGCGTGGCTGCTGGTGTGAGCTTCGTGTTCCGGGTGTAACCTGC
 TGCGCATCCTCCTGCCGTCCTGTCCGTTTATTTCCAGGGTGTGCAGAAGGATGCAGAGCAGGAGTC
 ACAGATGAGGGCAGAGATCCAGGGCATGAAGCAGGAGCTGTCCACCGTCAACATGATGGACGAGTTTGC
 AGATACGCCAGGCTGAAAGGAAGATCAACAAGATGACGGACAAGCTCAAACGCACGTGAAGGCACGGA
 CAGCTCAGTTGGCAAGATCAAGTGGTTTATAAGTGTGCTTTCTACATACTACAAGCAGCGCTGATGAT
 CTCGCTCATTGGAAATTACTCCGTCCTGTGGCTGTGGTGCCAGCAAGTGGATCACCCCACTAGAC
 CGCCTGGTAGCCTTTCCTACTCGAGTAGCAGGTGGAATTGGAATCACCTGTTGGATTTAGTCTGTAACA
 AAGTTGTGGCTATCGTTCTCCACCTTTTCAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201509 representing NM_207301
 Red=Cloning site Green=Tags(s)

MSASETDRWALLVLSFVFGCNLLRILLPSLSSFISRVLQKDAEQESQMRAEIQGMKQELSTVNMMEFA
 RYARLERKINKMTDKLKHVKARTAKLAKIKWFISVAFYILQAALMISLIWKYYSVPVAVVPSKWITPLD
 RLVAFPTRVAGGIGITCWILVCNKVVAIVLHPFS

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online >](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_207301

ORF Size: 522 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_207301.3</u>
RefSeq Size:	2483 bp
RefSeq ORF:	525 bp
Locus ID:	71446
UniProt ID:	<u>Q8K0D7</u>
Cytogenetics:	16 C4
Gene Summary:	Receptor for ASNA1/TRC40-mediated insertion of tail-anchored (TA) proteins into the ER membrane.[UniProtKB/Swiss-Prot Function]