

Product datasheet for **MG201530**

Reg3b (NM_011036) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Reg3b (NM_011036) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Reg3b
Synonyms: HIP; Pap; PAP1; REG-III
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201530 representing NM_011036
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCTGCCTCCAACAGCCTGCTCCGTCATGTCCTGGATGCTGCTCTCCTGCCTGATGCTCTTATCTCAGG
 TTCAAGGTGAAGACTCCCTGAAGAATATACCCTCCGCACGCATTAGTTGCCCAAGGGCTCCAGGCTTA
 TGGCTCCTACTGCTATGCCTTGTTCAGATACCACAGACCTGGTTTGATGCAGAAGTGGCTGCCAAAAG
 AGGCCTGGAGGACACCTCGTATCTGTGCTCAATAGCGCTGAGGCTTCATTCTTGCCTCCATGGTGAAGA
 GAACAGGAAACAGCTACCAATACACTTGGATTGGGCTCCATGACCCCACTCTGGGTGCAGAACCCAATGG
 AGGTGGATGGGAATGGAGTAACAATGACGTGATGAATTACTTTAACTGGGAGAGGAACCCATCTACTGCC
 TTAGACCGTGCTTCTGTGGCAGCTTGTCAAGAGCTTCTGGATTCTAAAATGGAGAGATATGACATGTG
 AGGTGAAGTTGCCCTATGTCTGCAAATTTACTGGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201530 representing NM_011036
 Red=Cloning site Green=Tags(s)

MLPPTACSVMSWMLLSCLMLLSQVQGEDSLKNIPSARISCPKGSQAYGSYCYALFQIPQTFDAELACQK
 RPGGHLVSVLNSAEASFLSSMVKRTGNSYQYTWIGLHDPTLGAEPNGGGWEWSNNDVMNYFNWERNPSTA
 LDRAFCSLSRASGFLKWRDMTCEVKLPYVCKFTG

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011036.1](#), [NP_035166.1](#)

RefSeq Size: 759 bp

RefSeq ORF: 528 bp

Locus ID: 18489

UniProt ID: [P35230](#)

Cytogenetics: 6 C3

Gene Summary: Bactericidal C-type lectin which acts against several intestinal Gram-positive and Gram-negative bacteria. Lacks antibacterial activity against S.typhimurium. May play a role in protection against infection with S.enteritidis by inhibiting its translocation from the gut lumen into intestinal tissues and further extraintestinal tissues.[UniProtKB/Swiss-Prot Function]