

Product datasheet for **MG201532**

Klra21 (BC064073) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Klra21 (BC064073) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Klra21
Synonyms: Klra-ps1; Ly49H; Ly49u
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201532 representing BC064073
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGCAAGTGACATCAAGTTAAAGGAAGAAATGTCGACAACCTGAGTCTATAGATTGTAGGCCAAGCAATG
AACTTCTGGAATACATCAAAAGAGAACAGGAGAGATGGAACAGTGAAAGCAAGACGGTTTTAGATTCCTC
ACGTGACACAGGCAGAGGTGTTAAATACTGGTTCTGCTATGGTACTAAATGTTATTATTTTCATCATGAAC
AAAACATCATGGAGTGGATGTAAAGCGAACTGCCAGCATTACAGCTTTCCATTGTGAAGATAGAAGATG
AAGATGAACTGAAATTCCTTCAACGCCATGTTATTCCAGAGAGTTACTGGATTGGATTGTCATATGATAA
GAAAAAAGGAATGGGCATGGATTGACAATGGCCCATCTAAACTTGACATGAAAAGGAAAATGAACTTT
ACGTCTAGAGGATGTGTATTTTATCTAAAGCAAGAATAGAAGATACTGACTGTAATACTCCCTACTACT
GTATTTGTGGGAAGAACTGGATAAATTCCTGAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201532 representing BC064073
Red=Cloning site Green=Tags(s)

MQSDIKLKEEMSTTESIDCRPSNELLEYIKREQERWNSSESKTVLDSSRDTGRGVKYWFCYGTKCYFFIMN
 KTTWSGCKANCQHYSFPIVKIEDEDELKFLQRHVIPESYIWGLSYDKKKKEWAWIDNGPSKLDMKRKMNF
 TSRGCVFLSKARIEDTDCNTPYYCICGKKLDKFPD

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




View online »

Cloning Scheme:


ACCN:	BC064073
ORF Size:	525 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC064073 , AAH64073
RefSeq Size:	1172 bp
RefSeq ORF:	527 bp
Locus ID:	93968

Cytogenetics: 6 F3

Product images:



Circular map for MG201532