

## Product datasheet for **MG201587**

### **Aipl1 (BC028285) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Aipl1 (BC028285) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Aipl1  
**Synonyms:** A930007I01Rik, MGC25485  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG201587 representing BC028285  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGACGTCTCTACTCCTCAATGTGGAGGGTGTCAAGAAGACCATTCTGCATGGGGGACGGGAGAGC  
 TCCCCAGCTTCATCACTGGCTCTCGAGTGACCTTTCATTTCCGCACAATGAAGTGTGATGACGAGCGCAC  
 GGTGATCGATGACAGCAAGCAGGTGGGCCAGCCATGAGCATCATCATCGGCAACATGTTCAAGCTGGAG  
 GTGTGGGAGACGCTGCTGACCTCCATGCGGCTGGGCGAGGTGGCTGAGTTCTGGTGGACACCATTACA  
 CAGGGGTCTACCCTATGTTGTCCCGAGTCTGCGGAGGTGGCTGAGGGCAAGGACCCACAAGCTGGCA  
 TGTGCACACGTGCGGGTTGGCCAACATGTTTGCATACCACACGCTGGGCTACGAGGACCTGGATGAGCTG  
 CAGAAAGAGCCACAGCCTCTTGTCTTCTGATTGAACTGTTGCAGGATTTTTATTGATTGGTCATCTGC  
 TTCCATTTCTGCAAGATTTAGAGTTGCGGCTCTTCTTGGTCCCTG

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >MG201587 representing BC028285  
 Red=Cloning site Green=Tags(s)

MDVSLLLNVEGVKKTILHGGTGELPSFITGSRVTFHFRTMKDDERTVIDDSKQVGPMSIIIGNMFKLE  
 VWETLLTSMRLGEVAEFWCDTIHTGVYPMLSRSLRQVAEGKDPTSHVHTCGLANMFAYHTLGYEDLDEL  
 QKEPQPLVFLIELLQDFYSIGHLLPFLQDLELRPLPWSL

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI




[View online >](#)

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** BC028285

**ORF Size:** 537 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC028285.1</a></u>
<b>RefSeq Size:</b>	2411 bp
<b>RefSeq ORF:</b>	539 bp
<b>Locus ID:</b>	114230
<b>Cytogenetics:</b>	11 B4
<b>Gene Summary:</b>	May be important in protein trafficking and/or protein folding and stabilization. [UniProtKB/Swiss-Prot Function]