

Product datasheet for **MG201601**

Aprt (NM_009698) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Aprt (NM_009698) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Aprt
Synonyms: C85684
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201601 representing NM_009698
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCGGAACCTGAGTTGAACTGGTGGCGCGGCATCCGCAGCTTCCCGACTTCCAATCCCGGGCG
 TGCTGTTCAAGGATATCTCGCCCTCTTGAAAGACCCGGACTCCTCCGAGCTTCCATCCGCCTTTGGC
 CAGTCACTGAAGTCCACGCACAGCGCAAGATCGACTACATCGCAGGTCTAGACTCCAGGGCTTCTGT
 TTTGGCCCTTCCCTAGCTCAGGAGCTGGCGTGGGCTGTGTGCTCATCCGAAACAGGGGAAGCTGCCGG
 GCCCACTGTGTCAGCCTCCTATTCTCTGGAGTATGGAAGGCTGAGCTGAAATCCAGAAAGATGCCTT
 GGAACCCGGGCAGAGAGTGGTCATTGTGGATGACCTCCTGGCCACAGGAGGAACCATGTTTGCGGCCTGT
 GACCTGCTGCACAGCTCCGGGCTGAAGTGGTGGAGTGTGTGAGCCTGGTGGAGCTGACCTCGTGAAGG
 GCAGGGAGAGGCTAGGACCTATACCATTCTCTCTCCTCCAGTATGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201601 representing NM_009698
 Red=Cloning site Green=Tags(s)

MSEPELKLVARIRSFDPFIPGVLFDRDISPLKDPDSFRASIRLLASHLKSTHSGKIDYIAGLDSRGFL
 FGPSLAQELGVGCVLIRKQKLPGPTVSASYSLEYGKAELEIQKDALEPGQRVVI VDDLLATGGTMAAC
 DLLHQLRAEVVECVSLVELTSLKGRERLGPFFSLLQYD

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_009698

ORF Size: 540 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_009698.2</u> , <u>NP_033828.2</u>
RefSeq Size:	866 bp
RefSeq ORF:	543 bp
Locus ID:	11821
UniProt ID:	<u>P08030</u>
Cytogenetics:	8 71.91 cM
Gene Summary:	Catalyzes a salvage reaction resulting in the formation of AMP, that is energetically less costly than de novo synthesis.[UniProtKB/Swiss-Prot Function]