

## Product datasheet for **MG201717**

### **Bmf (BC079650) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Bmf (BC079650) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Bmf  
**Synonyms:** AW260063; AW411872  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG201717 representing BC079650  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGAGCCACCTCAGTGTGTGGAGGAGCTAGAAGATGATGTGTTCCAGTCAGAGGATGGGGAGCCAGGGA  
 CACAGCCTGGGGCTTGCTCTGCTGACCTGTTTGCCAGAGCCAGCTGGACTGTCCCCTCAGTCGACT  
 CCAGCTCTTCCCTCACCCATTGCTGTGGTCCCGGACTCCGGCCATAAGCCAGGAAGACAAGGCCACT  
 CAGACCCTCAGTCCAGCTTCCCAAGCCAGGGTGTATGCTGCCTTGTGGGTGACAGAGGAACCCAGA  
 GACTCTTTTACGGCAACGCTGGCTACAGGCTTCTCTCCCTGCCAGTTTCCCTGCAGGCTCACCCCTTGG  
 GGAGCAGCCCTGAAGGACAGTTCCTTACGACCCGAGCAGAGGTGCAGATCGCCAGAAAGCTTCAGTGT  
 ATTGCAGACCAGTTCATCGGCTTCATACGCAACAACACCAGCAGAACCAGACCGTGCCTGGTGGCAGG  
 TCTTCTTTTCTTCAAACCTCGCCCTGAACAGACAAGAAAACAGGGAAGGGTGGGGCCCTGG

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >MG201717 representing BC079650  
 Red=Cloning site Green=Tags(s)


MEPPQCVEELEDVDFQSEDEPGTQPGLLSADLFAQSQDCPLSRLQLFPLTHCCGPGLRPISQEDKAT  
 QTLSPASPSQGVMLPCGVTEEPQRLFYGNAGYRLPLPASFPAGSPLGEQPPEGQFLQHRAEVQIARKLQC  
 IADQFHLHTQQHQNRDRAWVFLFLQNLALNRQENREGVGPW

**TRTRPLE** - GFP Tag - V


**Restriction Sites:** Sgfl-MluI



Cloning Scheme:



Plasmid Map:



ACCN: BC079650

ORF Size: 555 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC079650.1</a></u>
<b>RefSeq Size:</b>	4667 bp
<b>RefSeq ORF:</b>	557 bp
<b>Locus ID:</b>	171543
<b>Cytogenetics:</b>	2 E5
<b>Gene Summary:</b>	May play a role in apoptosis.[UniProtKB/Swiss-Prot Function]