

## Product datasheet for **MG201753**

### **Bcl2l12 (BC039807) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Bcl2l12 (BC039807) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Bcl2l12  
**Synonyms:** Bcl-L12, Bcl2-L12  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG201753 representing BC039807  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCCGGCTCTGAAGAGTTGGGTCTCCGAGAGGACACACTGAAGTTCTGACGGCCTTCCTTAAGCGCG  
 GTGAAGTCGCTGGGTCCCCTGTCCCACTCCACCCAGGAGTCTGCCAAGAGGAGACAACAGATTTCT  
 GAGCCGCCTCGAAGATGCCTCCCTGTCCCTGGGGAGAGGAGCTCCCCACGGAGTCTCTAGGCC  
 CACTTCTTCCCCTCGCCCTGCTATGGTTCGAGCCTGGCCAGCCTTCCGAGTCTATGCCCTGG  
 TTGCCAGCGGCTGGAACAATTGGTGCAGGAGCAACTGAAGTCCCGCCAGTTCAGAGTCCAGGGACC  
 CCCACCCACGGAGAAGGAAGCCCTGCTAAGGCGGCTGGTAGCCTTGTGGAGGAAGAGGCAGAAGTCATC  
 AACCAAAAGGTGACCAGCCTCTACTCCACGCTATGGAGAGGACAGTTGTGGAATAGGGGGCAGACGCTGG  
 CGCTGCTTTGTTTGTATGCTGGGTGTGGTTCAGGGCTCCAGTAAACATCAGAGCATCTTACAAGAGCC  
 T

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >MG201753 representing BC039807  
 Red=Cloning site Green=Tags(s)

MAGSEELGLREDTLKVLTAFLKRGEVAGSPVPTPPRSPAQEETDFLSRLRRCLPCPLGRGAPPTSSRP  
 HFLPLRPCYGSEPGPASSEFYALVAQRLEQLVQEQLKSPSSSEFQGGPPTEKEALLRRLVALLEEEAEVI  
 NQKVTSLSYTLWRGQLWNRGQTLALLCFDAGCGSRAPVNIRASSQEP

**TRTRPLE** - GFP Tag - V


**Chromatograms:** [https://cdn.origene.com/chromatograms/ja1755\\_a04.zip](https://cdn.origene.com/chromatograms/ja1755_a04.zip)



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC039807

ORF Size: 561 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC039807](#), [AAH39807](#)


RefSeq Size: 1620 bp

RefSeq ORF: 563 bp

Locus ID: 75736

Cytogenetics: 7 B3

Product images:



Circular map for MG201753