

## Product datasheet for **MG201798**

### Stmn4 (NM\_019675) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Stmn4 (NM\_019675) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Stmn4  
**Synonyms:** Lagl; R75359; RB3  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG201798 representing NM\_019675  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGACTCTTGCGCCTATAAGGAGAAGTGAAGGAACTCCCCTAGTGTCTCTGTTCTGCTCTTGCTTCT  
 TGTCTGATCCCCTGAATAAATCATCCTACAAATATGAAGCAGATACAGTAGACCTGAACTGGTGCAT  
 CTCTGATATGGAAGTCATCGAGCTGAACAAGTGCACCTCGGGCCAGTCTTTGAAGTCATCTTGAAGCCA  
 CCTTCTTTGACGGGTGCTGAGTTAATGCCTCCCTCCAAGCGCGAGACCCATCATTAGAAGAGA  
 TTCAGAAGAAGCTAGAAGCAGCAGAGGAGCGAAGGAAGTACCAGGAAGCTGAACTCCTGAAACACCTCGC  
 GGAGAAACGAGAGCACGAGCGAGAGGTGATCCAGAAAGCTATCGAAGAAAACAATAACTTCATCAAGATG  
 GCAAAAGAGAAGCTGGCCAGAAGATGGAGTCCAATAAGGAGAACCAGGGAGGCCATCTGGCTGCCATGT  
 TGGAGCGGCTGCAAGAGAAGGACAAGCACGCAGAGGAGGTGCGGAAAAACAAGGAGCTGAAGGAAGAGGC  
 CTCCAGG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG201798 representing NM\_019675  
 Red=Cloning site Green=Tags(s)

MTLAAAYKEKMKELPLVSLFCSCFLSDPLNKSSYKYEADTVDLNWCVISDMEVIELNKCTSGQSFEVILKP  
 PSFDGVPEFNASLPRRRDPSLEEIQKKLEAAEERRKYQEAELLKHLAEKREHEREVIQKAIENNFIK  
 AKEKLAQKMSNKENREAHLAAMLRLQEKDKHAEVVRKNKELKEEASR

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_019675

**ORF Size:** 567 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_019675.3](#), [NP\\_062649.1](#)

**RefSeq Size:** 1249 bp

**RefSeq ORF:** 570 bp


**Locus ID:** 56471

**UniProt ID:** [P63042](#)

**Cytogenetics:** 14 D1

Gene Summary: Exhibits microtubule-destabilizing activity.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MG201798