

Product datasheet for **MG201853**

Sec11c (NM_025468) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sec11c (NM_025468) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Sec11c
Synonyms: 1810029G24Rik; Sec1113
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201853 representing NM_025468
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTTCGTGCGGGCGCCGTGGGACGCATCTCCCACGTCAGCCTGGACATCTTTGGGGACCTGAGGA
 AGATGAACAAACGACAGCTCTACTACCAGGTTTTAACTTTGCCATGATCGTGTCTTCTGCGCTCATGAT
 CTGGAAAGGCCTGATTGTTCTCACGGGCAGCGAGAGTCCCATCGTGGTGGTACTCAGTGGCAGTATGGAG
 CCGGCCTCCACAGAGGAGATCTGCTGTCTCACGAATTTCCGGGAGGACCCCATCAGAGCTGGTGAAA
 TAGTTGTTTTTAAGGTTGAAGGAAGAGACATTCGATAGTTCACAGAGTAATCAAGGTTTCAAAAAAGA
 TAATGGTGACATCAAGTTTCTGACTAAAGGAGATAATAATGAAGTCGATGATAGAGGCTTGACAAAGAA
 GGCCGAACTGGCTGGAGAAGAAGGACGTGGTAGGAAGAGCTCGAGGGTCTTACCTTATGTTGGCATGG
 TCACCATCATAATGAATGACTATCCAAAATTCAAGTATGCTCTTTGGCTGTGATGGGTGCATATGTGTT
 ACTAAAACGTGAATCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201853 representing NM_025468
 Red=Cloning site Green=Tags(s)

MVRAGAVGTHLPTSSLDIFGDLRKMNRQLYYQVLFAMIVSSALMIWKGLIVLTGSESPIVVVLSGSME
 PAFHRGDLLFLTNFREDPIRAGEIVVFKVEGRDIPVHRVIKVHEKDNGDIKFLTKGDNNEVDDRGLYKE
 GQNWLEKKDVVGRARGFLPYVGMVTIIMNDYPKFKYALLAVMGAYVLLKRES

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_025468

ORF Size: 576 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_025468.2](#)

RefSeq Size: 771 bp

RefSeq ORF: 579 bp

Locus ID: 66286

UniProt ID: [Q9D8V7](#)

Cytogenetics: 18 E1

Gene Summary: Component of the microsomal signal peptidase complex which removes signal peptides from nascent proteins as they are translocated into the lumen of the endoplasmic reticulum.
[UniProtKB/Swiss-Prot Function]