

Product datasheet for **MG201864**

Csrp2 (NM_007792) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Csrp2 (NM_007792) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Csrp2
Synonyms: AW551867; Crp2; SmLim
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201864 representing NM_007792
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCTGTCTGGGGCGGTGAAATAAGTGC GGCCCTGCGGGAGAACCGTGTACCACGGGAAGAGGTGC
 AGTGGGATGGGCGGAGCTTCCATCGCTGCTTCTGTGCATGGTTGCAGGAAAAATTTAGACAGCAC
 AACAGTGGCGATTCATGATGAAGAGATCTACTGCAATCTGCTACGGAAAGAAGTATGGACCAAAAGGC
 TATGGTTATGGCAGGGCGCTGGCACGCTCAACATGGACCGGGTGGAGACTGGGCATCAAGCCAGAGA
 GTGCTCAACCTCACAGGCTACGACAAATCCAAACACTTCTAAATTTGCCAGAAATATGGAGGAGCTGA
 GAAGTGTCCAGGTGTGGGGATTCCGTGTATGCTGCGGAGAAGATCATTGGAGCTGGGAAGCCCTGGCAC
 AAAAAGTGTTCGGTGTGCAAGTGTGGGAAGAGTCTGGAGTCTACAACTCTGACTGAGAAAGAAGGCG
 AAATCTACTGTAAAGGTGCTACGCAAAGAAGCTTTGGGCCAAGGGATTTGGCTATGGTCAAGGGGCAGG
 GGCCCTTGTTTCATGCTCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201864 representing NM_007792
 Red=Cloning site Green=Tags(s)

MPVWGGGNKCGACGRTVYHAEVQCDGRSFHRCCFLCMVCRKNLDSTTVAIHDEEYCKSCYKGYGPKG
 YGYGQAGTLNMDRGERLGIKPESAPHRPTTNPNTSKFAQKYGGAEKCSRCGDSVYAAEKIIGAGKPWH
 KNCFRCAKCGKSLESTLTKKEGEIYCKGKYAKNFGPKGFGYGGAGALVHAQ

TRTRPLE - GFP Tag - V


Restriction Sites: Sgfl-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM_007792

ORF Size: 579 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_007792.4](#)

RefSeq Size: 995 bp

RefSeq ORF: 582 bp

Locus ID: 13008

UniProt ID: [P97314](#)

Cytogenetics: 10 D1

Gene Summary: Drastically down-regulated in response to PDGF-BB or cell injury, that promote smooth muscle cell proliferation and dedifferentiation. Seems to play a role in the development of the embryonic vascular system (By similarity).[UniProtKB/Swiss-Prot Function]