

Product datasheet for **MG201928**

Rer1 (NM_026395) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rer1 (NM_026395) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Rer1
Synonyms: 1110060F11Rik; 5830454N22Rik; AU043380
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201928 representing NM_026395
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCTGAAGGGGACAGTGTGGAGACTCTGTCCATGGGAAGCCCTCGGTGGTGTACAGATTTTTCTCAC
 GGCTTGGACAGATTTATCAGTCTGGCTGGACAAGTCTACCCCTACACAGCCGTCGGTGGGTTGTGAC
 ACTGGCCTGAGTTTTGTCTACATGATTAGAGTTACCTGTTACAAGTTGGTACATTGTGACCTACGCC
 TTGGGAATTTACCACCTAAACCTTTTCATTGCGTTCCTTTCTCCAAAGTGGATCCTTCTTGATGAAG
 ACTCAGATGACGGCCCTCATTACCAACCAACAGAATGAGGAATTCGTCCTTCATTGAAGGCTTCC
 AGAATTCAAATTTGGCAGCAGCTACAAAAGGCATTCTGTGGCTATGATCTGTACCTTCTCGAGGCT
 TTTAACGTCCCAGTGTCTGGCCGATCCTGGTGTGACTTCATCATGCTTTTCTGTATCACAATGAAGA
 GGCAAATCAAGCACATGATTAATACCGGTACATACCATTACACACGGGAAGAGGAGATACAAAGGGAA
 GGAGGACGTGGCAAGACATTTGCTAGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201928 representing NM_026395
 Red=Cloning site Green=Tags(s)

MSEGDSVGDVSHGKPSVYRFFSRLGQIQSWLDKSTPYTAVRWVTLGLSFVYMIRVYLLQGWIVTYA
 LGIYHLNLFIAFLSPKVDPSLMEDSDGPSLPTKQNEEFRPFIRRLPEFKFWHAATKGILVAMICTFFEA
 FNVVPVFWPILVMYFIMLFCITMKRQIKHMIKYRYPFTHGKRRYKGEDVVKTFAS


TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_026395

ORF Size: 588 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026395.2](#)

RefSeq Size: 1518 bp

RefSeq ORF: 591 bp

Locus ID: 67830

UniProt ID: [Q9CQU3](#)

Cytogenetics: 4 E2

Gene Summary: Involved in the retrieval of endoplasmic reticulum membrane proteins from the early Golgi compartment.[UniProtKB/Swiss-Prot Function]