

Product datasheet for **MG202026**

1700001C02Rik (NM_029285) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700001C02Rik (NM_029285) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: 1700001C02Rik
Synonyms: 1700047H18Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG202026 representing NM_029285
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCTTCCGTAGCGCGGGCACACTAATGACCGAATTTAATGCCGCTTCGTGCCCCAGCTCTCATGC
 CCGGGTACAAAGGCCATGTTCCCGGTGTGGCTTCTCCTTTGGCTCTTCTATGGCAACACCACCTTCAA
 GTACTTCCAGGACCTCCGCAACACAGGCTTGGAGAAGAGCTATGCGCTCCTCAGTGGGGTTGCTTTCCC
 ACCATCTTCTCCCCTAACTCCCCCTGGTCTGACTGATAACTCACAGAAGCTGGGATCGATGGCTGCACC
 AGCCACCTACAGCCGATCCAACCAGGACAACGGGCGCACCAGTACTACTACTTTTACAAGACAGT
 CCAGCAGCAGCGAAGTATTACCAAGACAAGACTGGCATTCTGTCCCGTGTCTTACTTTGACTTCTCT
 GTGAAGGAGCCGACCGGTACCCCTACCCACCGACCTGCCTCCTCTGAGTGAGAAGAGTAAAGTGAACC
 TACTAAGATTATCCCCTGACTACAAGAGGACCTACCAGACCTTCCCATCAGGGAAGAGAGTTTCTCTCA
 AGAGCGCCAAAGGCGGGACCTTTACTTCGAGTTCAGAGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG202026 representing NM_029285
 Red=Cloning site Green=Tags(s)

MAFRSAGTLMTEFNAAFVPPALMPGYKGVHVPVASSFGSSYGNTTFKYFQDLRNTGLEKSYALLSGGCFP
 TIFSPNSPLVLTDNSQNWDRLHQPTYSRNQDNGRTELTFYKTVQQQRKYQDKTGILSRVPYFVLP
 VKEPDRYPLPTDLPLSEKSKWNLLRLSPDYKRTYQTFPSGKRVSSQERQRDLVFEFRA

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_029285

ORF Size: 600 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_029285.1](#), [NP_083561.1](#)

RefSeq Size: 712 bp

RefSeq ORF: 603 bp

Locus ID: 75434

UniProt ID: [Q9DAS2](#)

Cytogenetics: 5 B1