

Product datasheet for **MG202035**

Commd7 (BC094216) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Commd7 (BC094216) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Commd7
Synonyms: 2310010I22Rik; AW060956; mU3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG202035 representing BC094216
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGCCCGCTGCACTGCACGCAGGACCCGGTGCCCGAGGCCGTGCGCGGAGACATGCAGCAGCTGAACC
 AGCTGGGCGCTCAGCAGTTCTCCGACCTGACAGAAGTGCTTTTCCATTTCCCTAACTGAGCCCAAAGAGGT
 GGAGCGATTTCTGGCTCAGCTCTCTGAATTTGCCACCAGCAATCAGATCAGTCTTGGCCCCCTCAAAGC
 ATCATGAAAAGCCTCCTACTGGTCCAAATGGTGCCTGAAGAAAGCCCTCACTGCCGAGCAGGTGCGGA
 CTGATCTCCAGACGCTAGGTCTTAGCGAGGAGAAAGCCACTTACTTTTCTGAAAAGTGGAAGCAGAACGC
 CTCGACCCTTGACAGTGGGCCATGGGTGACACCTGATGGTTAACCCAGCTGATTGATATGGAGTGGAGG
 TTTGGAGTGACGTCTGGGAGCAGTGAATTAGAGAAAGTGGGAAGCATATTTTACAGTTAAAGTTGGTGG
 TTAAGAAAGGAAAACAACTGAAAATTTGTATATGGAACCTAACCTTGCCCGAGTTCTACAGCTTCTCGCA
 CGAGATGGAGAGAGTCCGAGCGAGCATGGAGTGCCTCAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG202035 representing BC094216
 Red=Cloning site Green=Tags(s)

MGRLHCTQDPVPEAVRGDMQQLNQLGAQQFSDLTEVLFHFLTPEPKEVERFLAQLSEFATSNQISLGPLKS
 IMKSLLLVPNGALKKGLTAEQVRTDLQTLGLSEEKATYFSEKWKQNASTLAQWAMGQTLMVNQLIDMEWR
 FGVTSGSSELEKVGSIIFLQLKLVVKKGKQENLYMELTLPQFYFSLHEMERVRASMECLS

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI



Cloning Scheme:



Plasmid Map:



ACCN: BC094216

ORF Size: 602 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC094216](#), [AAH94216](#)

RefSeq Size: 1287 bp

RefSeq ORF: 602 bp

Locus ID: 99311

Cytogenetics: 2 H1

Gene Summary: May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes. Associates with the NF-kappa-B complex and suppresses its transcriptional activity.[UniProtKB/Swiss-Prot Function]