

Product datasheet for **MG202186**

Tspan15 (BC003872) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tspan15 (BC003872) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Tspan15
Synonyms: 1110036D12Rik; 2700063A19Rik; AW048364; Tm4sf15
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG202186 representing BC003872
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCTCGCGGGGACTCGGAGCAGGTGCGCTACTGCGCGCGTTCTCCTATCTTTGGCTCAAGTTCTCTC
TCATCATCTACTCCACCGTGTCTGGCTGATTGGGGCCCTGGTCTGTGGTGGGATCTACGCAGAGGC
AGAGCGGCAGAAATACAAAACCCTGGAGAGTGCCTTCTGGCCCCGCCATTATCCTCATCTCTGGGG
GTGGTCATGTTTCATCGTCTCCTTCATCGGGGTGCTGGCTCCCTCGGGACAACCTGTGCCTTCTGCAGT
CGTTTATGTATATCCTGGGGATCTGCCTGGTCATGGAGCTGATTGGTGGCATAAGTGGCTTTGATCTCCG
GAACCAGACTATTGACTTCTGAACGACAACATTCGGAGAGGAATCGAGAATTACTACGATGATCTGGAC
TTCAAGAACATCATGGACTTTGTTCAAGAAGTTCAAGTGTGTGGCGGGGAGGACTACAGAGACTGGA
GCAAAAACAGTACCACGACTGCAGCGCCCCGGGCCCTGGCTGCGGGGTGCCCTACACCTGCTGCAT
CAGGAACACGATGTTGTCAACACCATGTGTGGCTACAAAACAATCGACAAGGAGCGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG202186 representing BC003872
Red=Cloning site Green=Tags(s)

MPRGDSEQVRYCARFSYLWLFSLIIYSTVFWLIGGLVLSVGIYAEAERQKYKTLES AFLAPAILILLG
VVMFIVSFIGVLAASLRDNLCLLQSFMYILGICLVMELIGGIVALIFRNQITDFLNDNIRRG IENYYDDL
FKNIMDFVQKKFKCCGGEDYRDWSKNQYHDCSAPGLACGVPYTC CIRNTMLSTPCVATKQSTRSA

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: BC003872

ORF Size: 620 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC003872](#), [AAH03872](#)

RefSeq Size: 1661 bp

RefSeq ORF: 620 bp

Locus ID: 70423

Cytogenetics: 10 B4

Gene Summary: Regulates maturation and trafficking of the transmembrane metalloprotease ADAM10 (PubMed:23035126, PubMed:26668317). Promotes ADAM10-mediated cleavage of (CDH2) (PubMed:26668317). Negatively regulates ligand-induced Notch activity probably by regulating ADAM10 activity (By similarity).[UniProtKB/Swiss-Prot Function]