

## Product datasheet for **MG202237**

### Mpzl1 (NM\_001001880) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Mpzl1 (NM\_001001880) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Mpzl1  
**Synonyms:** 1110007A10Rik; PZR  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG202237 representing NM\_001001880  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCAGAGGCCGTGGAGCCGTGGCACTGATTGCGGCCCGGCCCGCCGGCGCTGGCTGTGGTCCGGTGC  
 TAGCCCGATGCTCGGGCTGTTGACAGCTAGAATATCAGCCTGGAGGTCCATACTCCTAAAGAAATCTT  
 TGTGGTGAATGGGACACAAGGGAAGCTGACTTGCACATTCGATTCTCCTAACACAACCTGGATGGTTGACC  
 ACGGTCTCCTGGAGCTTCCAGCCAGACGGCACCAGACGCGCGTGTGTTTTTCCACTACTACAAGGAC  
 AGGTGTACATTGGGGATTATCCACCATTTAAAGACCGAGTCACCTGGGCTGGGGACCTTGACAAGAAAGA  
 TGCATCAATCAATATAGAAAATATTCAGGCTGTTCACAATGGCACCTATATCTGGCATGTCAAAAATCCT  
 CCTGACATTGTTGTGCGCCTGGACACATTAGGCTCCACGTGGTGGAAATAGACAACCTGTTGGTGTTC  
 TGGTTTGGGTGGTGGGCACTGTCACTGCTGTGGTCTTGGCCTCACTCTGCTCATCAGCTTGGTCTCT  
 GGTGCTCTCTACAGAAGGAAACATTCGAAGCGGGATTATACCGGGGCCAGTCATTTACGCACAGT

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >MG202237 representing NM\_001001880  
 Red=Cloning site Green=Tags(s)

MAEAVGAVALIAAPARRRWLWSVLAAMLGLLTARISALEVHTPKEIFVVNGTQGKLTCTFDSPNTTGWLT  
 TVSWSFQPDGDSAVSFFHYSQGVYIGDYPPFKDRVTWAGDLDDKSDASINIENIQAVHNGTYICDVKNP  
 PDIVVRPGHIRLHVVEIDNLLVFLVWVVVGTAVVLGLTLLISLVLVLYRRKHSKRDTGAQSFTHS

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001001880

ORF Size: 627 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001001880.2](#)

**RefSeq Size:** 2205 bp

**RefSeq ORF:** 630 bp

**Locus ID:** 68481

**UniProt ID:** [Q3TEW6](#)

**Cytogenetics:** 1 H2.3

**Gene Summary:** Cell surface receptor, which is involved in signal transduction processes. Recruits PTPN11/SHP-2 to the cell membrane and is a putative substrate of PTPN11/SHP-2. Is a major receptor for concanavalin-A (ConA) and is involved in cellular signaling induced by ConA, which probably includes Src family tyrosine-protein kinases. Isoform 2 seems to have a dominant negative role; it blocks tyrosine phosphorylation of MPZL1 induced by ConA. Isoform 1, but not isoform 2, may be involved in regulation of integrin-mediated cell motility (By similarity).[UniProtKB/Swiss-Prot Function]