

Product datasheet for **MG202295**

Timp3 (NM_011595) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Timp3 (NM_011595) Mouse Tagged ORF Clone
 Tag: TurboGFP
 Symbol: Timp3
 Synonyms: Timp-3
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-AC-GFP (PS100010)
 E. coli Selection: Ampicillin (100 ug/mL)
 ORF Nucleotide Sequence: >MG202295 representing NM_011595
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGACTCCCTGGCTTGGGCTTGTCTGTGCTCCTGAGCTGCTGGAGCCTTGGGCACTGGGGCGCGGAAGCGT
 GCACATGCTCTCCAGCCACAGCCAGGATGCCTTCTGCAACTCCGACATCGTGATCCGGGCCAAAGTGGT
 GGGAAAGAAGCTGGTGACGGAGGGGCCCTTGGCACTCTGGTCTACACTATTAAGCAGATGAAGATGTAC
 CGAGGCTTCAGTAAGATGCCCCACGTGCAGTACATTCACACGGAAGCCTCTGAAAGTCTTTGTGGCCTCA
 AGCTAGAAGTCAACAAATACCAGTACCTGCTGACAGGGCGCGTGTATGAAGGCAAGATGTACACAGGACT
 GTGCAACTTTGTGGAGAGGTGGGACCACCTCACACTGTCCAGCGCAAGGGCCTCAATTACCGCTACCAC
 CTGGGTTGCAATTGCAAGATCAATTCCTGCTACTACTTGCCTTGTGTTTGTGACCTCCAAGAATGAGTGTC
 TCTGGACCGACATGCTCTCCAATTTGGGTACCCTGGCTATCAGTCCAAACACTACGCCTGCATCCGGCA
 GAAGGGTGGTACTGCAGCTGGTACCGAGGATGGGCTCCCCAGACAAGAGCATCAGCAACGCCACAGAC
 CCC

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG202295 representing NM_011595
 Red=Cloning site Green=Tags(s)

MTPWLGLVLLSCWSLGHWGAEACTCSPSHSQDAFCNSDIVIRAKVVGKLVTEGPFGLVYTIKQMKMY
 RGF SKMPHVQYIHTEASESLCGLKLEVNKYQYLLTGRVYEGKMYTGLCNFVERWDHLTLSQRKGLNYRYH
 LGCNCKINSCYLLPCFVT SKNECLWTDMLS NFGYPGYQSKHYACIRQGGYCSWYRGWAPPDKSISNATD
 P

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_011595


ORF Size: 633 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_011595.2</u> , <u>NP_035725.1</u>
RefSeq Size:	4681 bp
RefSeq ORF:	636 bp
Locus ID:	21859
UniProt ID:	<u>P39876</u>
Cytogenetics:	10 42.83 cM
Gene Summary:	Complexes with metalloproteinases (such as collagenases) and irreversibly inactivates them by binding to their catalytic zinc cofactor. May form part of a tissue-specific acute response to remodeling stimuli.[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for MG202295