

## Product datasheet for **MG202311**

### Rab15 (NM\_134050) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Rab15 (NM_134050) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Rab15
Synonyms:	2310012G06Rik; AI840042
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG202311 representing NM_134050 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGAAACAGTACGATGTGCTGTTCCGGCTACTGCTGATCGGGGACTCCGGGTTGGCAAGACATGCC  
TGCTGTGCCGCTTCACCGACAACGAGTTCCTCCTCGCATATCTCCACCATCGGTGTTGACTTTAAGAT  
GAAGACTATCGATGTAGACGGCATCAAAGTGAGAATACAGATTTGGGACACAGCAGGGCAGGAGAGGTAC  
CAGACTATCACAAGCAGTACTATCGGCGAGCCCAGGGAATATTTTAGTCTACGACATTAGCAGTGAGC  
GCTCCTATCAGCATATCATGAAGTGGGTCAGTGACGTGGATGAGTACGCTCCAGAAGGAGTCCAGAAGAT  
CCTAATTGGGAATAAGGCTGATGAAGAGCAGAAACGGCAGGTGGGAGAGAGCAGGGGCAGCAGCTGGCT  
AAGGAGTACGGCATGGACTTCTACGAAACAAGTGCCTGCACCAACCTTAATATTAAGAGTCTTCACTC  
GTCTGACGGAGCTGGTGTGCTGCAGGCCACAGGAAAGAGCTGGATGGTCTCCGAACACGTGCCAGCAACGA  
GCTCGCACTGGCCGAGCTGGAGGAGGACGAAGGCAAACTGAGGGCCAGCAAACCTTCAAAGACCTGC  
TGGTGC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG202311 representing NM\_134050  
Red=Cloning site Green=Tags(s)

MAKQYDVLFRLLLIGDSGVGKTCLLCRFTDNEFHSSHISTIGVDFKMKTI DV DGIK VRIQI WDTAGQERY  
 QTITKQYYRRAQGI FLVYDISSERSYQHIMKWVSDVDEYAPEGVQKILIGNKADEEQKRQVGREQQQLA  
 KEYGMDFYETSACTNLNIKESFTRL TELVLQ AHRKELDGLRTRASNELALAELEEDEGKPEGPANSSKTC  
 WC


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_134050

**ORF Size:** 636 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_134050.3</a>
<b>RefSeq Size:</b>	3139 bp
<b>RefSeq ORF:</b>	639 bp
<b>Locus ID:</b>	104886
<b>Cytogenetics:</b>	12 C3