

Product datasheet for **MG202374**

Fam167a (NM_177628) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fam167a (NM_177628) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Fam167a
Synonyms: A030013D21
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG202374 representing NM_177628
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTGTGCCCCAGATCCAAGTGAAGAAGTGAGTGAGAAAGACAGGCCAGCCGGGGCAGCCGTACCTC
CTGATGACCACCTCCTCAGTCTGAAGGCCCTCACTGAGAAGCTGAGGCTGGAGACCCGAAGACCTTCCTA
CCTGGAATGGCAGGCCAGGCTGGAAGAGCAGACATGGCCTTCCCAAGGCCAGCTGCCAGCAAGAGGCG
AGTCTGGAGCAGGGAGCATGTGGGGTGGAGAGCCTTTGATGCCTCTAAAGGAGCCCGAGACCTTCTCC
CACCTTCCGCGAGTGCCGCGAGGGGTGACAGGCCACTGACCACGGCAAGCTGGAAGGCTTCCAGAGTAT
CGATGAAGCTATAGCCTGGCTCAGAAAGGAAGTGGCAGAGATGCGGCTCCAGGACCAGCAGCTGGCCAGG
CAGCTCATGCGCCTCCGAGGAGACATCAACAAGTTGAAGATCGAGCAGACCTGCCGCTGCACCGGAGGA
TGCTCAACGATGCCGCTTTGAGCTGGAGGAGCGGGATGAGCTGTCCGACCTCTTTTGTGACTCCCCGCT
GGCATCTCTTTTCCAGCCTTTCCATGCCACTCAAGCTCATTGGCGTCACCAAAATGAACATCAACTCCAGA
AGGTTCTCTCTGTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG202374 representing NM_177628
 Red=Cloning site Green=Tags(s)

MSVPQIQVEEVSEKDRPAGAAVPPDDHLLSLKALTEKLRLRTRRPSYLEWQARLEEQTWPFPRPAAQEEA
 SLEQGACGGGEPLMPLKEPRDLLPPSASAGRDRPLTTGKLEGFQSIDEIAWLKELAEMLRQDQQLAR
 QLMRLRGDINKLKEIQTCLHRRMLNDAAFELEERDELSDLFCDSPASSFSLSMPLKLIQVTKMINSR
 RFSLC


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_177628

ORF Size: 645 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_177628.4 , NP_808296.2
RefSeq Size:	4033 bp
RefSeq ORF:	648 bp
Locus ID:	219148
UniProt ID:	Q6P1G6
Cytogenetics:	14 D1