

Product datasheet for **MG202394**

Kcnip1 (NM_027398) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Kcnip1 (NM_027398) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Kcnip1
Synonyms: KCHIP1; Kchip1.2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG202394 representing NM_027398
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGGCCGTCATGGGCACTTCTCCTCCCTGCAGACCAACAAAGGCGACCCTCTAAAGACAAGATTG
AGGATGAGCTAGAGATGACCATGGTTTGCACCGGCTGAGGGACTGGAGCAGCTTGAGGCACAGACGAA
CTTACCAAGAGAGAAGTCAAGTCTTGTACCGGGATTCAAAAACGAGTGCCCTAGCGGTGTGGTCAAT
GAAGAAACATTCAAGCAGATCTACGCTCAGTTTTTCCCTCACGGAGATGCCAGCACATATGCACATTACC
TCTTCAATGCCTTCGACACCACCCAGACAGGCTCTGTAAGTTTCGAGGACTTTGTGACTGCTCTGTCGAT
TTTACTGAGAGGGACAGTCCATGAAAACTAAGGTGGACGTTTAAATTTGTATGACATCAATAAAGACGGC
TACATAAACAAAGAGGAGATGATGGACATAGTCAAAGCCATCTATGACATGATGGGAAATACACCTATC
CTGTGCTCAAAGAGGACACTCCCAGGCAGCATGTGGATGTCTTCTCCAGAAAATGGATAAAAATAAAGA
TGGCATTGTAACGTTAGATGAATTTCTTGAATCATGTCAGGAGGATGACAACATCATGAGATCTCTACAG
CTGTTCCAAAATGTCATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG202394 representing NM_027398
Red=Cloning site Green=Tags(s)

MGAVMGT FSSLQTKQRRPSKDKIEDELEMTMVCHRPEGLEQLEAQTNFTKRELQVLYRGFKNECPSGVVN
 EETFKQIYAQFFPHGDASTYAHYLFNAFDTTQTGSVKFEDFVTALSILLRGTVHEKLRWTFNLYDINKDG
 YINKEEMMDIVKAIYDMMGKYTPVLKEDTPRQHV DVFVQKMDKNKDGIVTLDEFLESCQEDDNIMRSLQ
 LFQNV M


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_027398

ORF Size: 648 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_027398.4
RefSeq Size:	1783 bp
RefSeq ORF:	651 bp
Locus ID:	70357
UniProt ID:	Q9JJ57
Cytogenetics:	11 A4
Gene Summary:	Regulatory subunit of Kv4/D (Shal)-type voltage-gated rapidly inactivating A-type potassium channels. Regulates channel density, inactivation kinetics and rate of recovery from inactivation in a calcium-dependent and isoform-specific manner. Modulates KCND2/Kv4.2 currents (PubMed:14572458). In vitro, modulates KCND1/Kv4.1 currents (By similarity). Increases the presence of KCND2 at the cell surface.[UniProtKB/Swiss-Prot Function]