

## Product datasheet for **MG202407**

### 2700062C07Rik (NM\_026529) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Tag:** TurboGFP  
**Symbol:** 2700062C07Rik  
**Synonyms:** A1195775; C87515  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)

**ORF Nucleotide Sequence:** >MG202407 representing NM\_026529  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTACTATAGGGCGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGCGGCAGAAGTATTACATTGAAGCGGCTGCTCGGGCCTGGTGGGAAGCTGCCCGGCCAGGCTCGGT  
 ACCTCTATGGGCCTACAGTTCGACACACGAGGATAACAGCACTTCCAGGAAACATGTCCACACTGCTT  
 CCAATTGTTGGTCTGGATAACTCTAGAGTGCGCCTCAAACCTAAAGCCAAATGACACCAAAGATACAG  
 AAATCTTAAATCGGGAAGCAAGAAATGATACACTCAGTTTTAAAGAAGCAAACTGTTGAGAAAAGTATA  
 AAGACTCTACAAGTGTGCTGCTGATTACCTGTAGAACATGCAACAGAACAGTGAGACACCATGGTAAGAG  
 TAGAAGCTTCTGTGGGCTCTGAAGAGCAACGCCGCACTGCTGCCAACAAAGCCAGCCGAAGACGCCA  
 AAGAGAACGGCTCCAGGCTCTGCAAACCTCGGTCAGAGTACAAATGGTTCCAAAGGCAAGAGTCCCTCCT  
 TGACCATCAGAACACCTACATCTGGACAGTCAACACCCATTTGCTCCTCAAGGAATGGGAGCAAAAGAAA  
 GAAACACTTCTCTCAACTAAAAGCACTGCTTAGTCAGAGTGATCTGACAAGAACCCAAAGCTGGACTTC  
 AGGCACTTCTTATCTTCTCTA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG202407 representing NM\_026529  
 Red=Cloning site Green=Tags(s)

MRQKYIIEAAARGLVGSCPGQARYLLWAYSSTHEDNSTFQETCPHCFQLLVLDNSRVRLKPKAKLTPKIQ  
 KLLNREARNDLSFKEAKLLRKYKDSVLLITCRTCNRTVRHHGKSRFSLWALKSNAATAANKASPKTP  
 KRTAPGSANLQSTNGSKGKSPSLTIRTPSGQSTPICSSRNGSKRKKHFSQLKALLSQSASDKNPKLDF  
 RHFLSSL

TRTRPLE - GFP Tag - V



**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_026529

**ORF Size:** 651 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_026529.3](#)

RefSeq Size: 1097 bp  
 RefSeq ORF: 654 bp  
 Locus ID: 68046  
 UniProt ID: [Q5XFZ0](#)  
 Cytogenetics: 18 A2

**Product images:**



Circular map for MG202407