

Product datasheet for **MG202474**

Sly (BC049626) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sly (BC049626) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Sly
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG202474 representing BC049626
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTCTTAAGAAATTGAAGGTGATACCAAAGGAAGTTACTTATTACTTTTGGACTTTGATGATGAGG
ACGATGACATAAAAGTTTCAGAGGAGGCTCTTTCGGAAGTAAAGAGCCCAGCATTGATAAAAATGAGAA
TATATCGCCTCAAGCAGAAGCAGATGAAGATATGGGAGATGAAGTAGACAGTATGTTGGATAAATCTGAA
GTAAATAATCCAGCAATTGGTAAAGATGAAAATATATCGCCTCAAGTAAAAGGAGATGAAGACATGGGAC
ATGAAGTAGGCAGTATGTTGGATAAATCTGGAGATGACATTTATAAGACGCTTCACATAAAGAGAAAAATG
GATGGAAACTTATGTCAAAGAATCTTTCAAAGGCAGCAACCAGAAATTAGAAAGATTTTGCAAAACGAAC
GAACGAGAGAGGAAGAACATCAACAACAAATTTGTGAGCAGTATATAACTACATTTCAGAAGTCTGATA
TGGATGTACAGAAATTCATGAAGAAAAAGAAAAATCAGTGAATAGTTGTCAAAAAGAACAACAAGCATT
GAAACTGTCCAAATGTAGTCAGAACCCAGACCCTGGAAGCAGTTAAAGAAATGCATGAGAAGTCCATGGAG
GTTTTGATGAACCTGGGGACCAAGAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG202474 representing BC049626
Red=Cloning site Green=Tags(s)

MALKKLKVIPKEGYLLLLDFDDEDDDIKVSEEALSEVKSPAFDKNENISPQAEADEMDGDEVDSMLDKSE
VNNPAIGKDENISPVKGDEDMGHEVGSMGDKSGDDIYKTLHIKRWMETVYKESFKGSNQKLERFCKTN
ERERKNINNKFCQYITTFQKSDMDVQKFNEEKEKSVNSCQKEQQALKLSKCSQNTLEAVKEMHEKSME
VLMNLGTFKN

TRTRPLE - GFP Tag - V




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



EcoRI BamHI KpnI RBS Kozac Consensus SgfI AscI
 CTATAGGCGCGCGGAATTCGCTGACTGGATCCGGTACCGAGAGATCTGCCGCCGATCGCCGCGCGCAGATCT
 HindIII NheI RsrII MluI NotI XhoI GFP Tag
 CAAGCTTAAGTAGCTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC ---
 T R T R P L E M E S D - - -
 PmeI FseI
 --- GAA GAA AGA GTT TAA ACGCCGGCCGCGAGCT
 - - E E R V Stop

ACCN: BC049626

ORF Size: 659 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC049626](#), [AAH49626](#)


RefSeq Size: 946 bp

RefSeq ORF: 659 bp

Locus ID: 382301

Cytogenetics: Y C2

Product images:



Circular map for MG202474