

Product datasheet for **MG202478**

Rpain (NM_027186) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rpain (NM_027186) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Rpain
Synonyms:	2400006N03Rik; Rip
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG202478 representing NM_027186 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGAGTCCTCGGGTCTCCGCACCGCTTGTGTACAAGCAGGTGGGCTCGCCCCACTGGAAAGAAA
CTTTCAGGCAGGGATGTCTGGAGAGAATGAGAAACAGCAGGCACAGGCTCCTGAACAAATATCGCCAGGC
TGCAGGTAGCACGCCGGGACAGCCTCAGACAGACTTCTGTGCAAGAAGTAATGGAGGAAGAGTGGGCT
TCTTTCAGTCTGTGGAGAATTGTCCGGAGGCCTTGCTTCAGTTGGAATTGCCACTGGACCTAGCTGTGC
TGCAGGACATCGAGCAGGAGCTGTGTAATGAAGAAAAGTCCATCATAAGTGAGTACGAAGAGGACTTAGA
GTTTGATGAAAGTTGTCTCAGGAGAATGTTGGCTGAGTGGGAAGCAAACCTCCCTCATCTGTCTGTGTGT
ATAAAGTACAACCTGAGAATAATGAACAGTGTGGTCACGTGTCCATGTGGCCTGCACATCCCTGTTCACT
CAACAGACCTGACAGAGCAGAAGCTTCGAGCCTGTTTGGAAAGAAAATGTGAATGAGCACAGTGTACACTG
TCCCCACACCCCTGTGTTCTCAGTCACCGGTGGAACAGAAGAGAAGCCAGTCTTCTGATGAACTGCCTG
ACTTGTGACACCTGGGCTGTGATCCTC

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG202478 representing NM_027186
 Red=Cloning site Green=Tags(s)

MAESSGSPHRLLYKQVGSPhWKETFRQGCLERMNRNRHLLNKYRQAAGSTPGTASDRLLVQEVMEEEWA
 SLQSVENCPEALLQLELPLDLAVLQDIEQELCNEEKSIISEYEEDLEFDESCLRRMLAEWEANSLICPVC
 IKYNLRIMNSVVTCPGLHIPVHSTDLTEQKLRACLEENVNEHSVHCPTPVF SVTGGTEEKPSLLMNCL
 TCDTWAVIL


TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_027186

ORF Size: 657 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_027186.1 , NP_081462.1
RefSeq Size:	1093 bp
RefSeq ORF:	660 bp
Locus ID:	69723
UniProt ID:	Q9CWY9
Cytogenetics:	11 B4
Gene Summary:	Mediates the import of RPA complex into the nucleus, possibly via some interaction with importin beta. Sumoylation mediates the localization of RPA complex into the PML body of the nucleus, thereby participating in RPA function in DNA metabolism (By similarity). [UniProtKB/Swiss-Prot Function]