

Product datasheet for **MG202521**

Gsta3 (NM_001077353) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gsta3 (NM_001077353) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Gsta3
Synonyms:	Gst2-3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG202521 representing NM_001077353 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGGGAAGCCAGTCCTTCATTACTTTGATGGCAGGGGAAGAATGGAGCCTATCCGGTGGCTCTTGG
CTGCAGCTGGTGTGGAGTTTGAAGAAAAATTTCTGAAACTCGGGATGACCTGGCAAGGTTACGAAGTGA
TGGGAGTCTGATGTTCCAGCAAGTGCCCATGGTAGAGATCGACGGGATGAACTGGTGCAGACCAAAGCC
ATTCTCAACTACATTGCCTCAAATACAACCTCTATGGGAAGGACATGAAGGAGAGGCCATCATTGACA
TGTACACAGAAGGAGTGGCGGATCTGGAGATAATGATTCTCTATTACCCACATGCCCCCTGAGGAGAA
AGAGGCAAGCCTTGCCAAGATCAAGGAACAACCAGGAACCGTTACTTCCCTGCCTTTGAAAAGGTGTTG
AAGAGCCATGGACAAGATTATCTCGTTGGCAACAGGCTGAGCAGGGCTGATATTGCCCTGGTTGAACTCC
TCTACCATGTGGAAGAGCTGGACCCGGCGTGTGGACAACCTTCCCTCTCCTGAAAGCGCTGAGAAGCAG
AGTCAGCAACCTCCACAGTGAAGAAGTTTCTCAACCTGGCAGCCAGAGGAAGCCTTTTGATGACGCA
AAATGTGTTGAGTCAGCAAAGAAGATTTTCAGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG202521 representing NM_001077353
 Red=Cloning site Green=Tags(s)

MAGKPVLYHYFDGRGRMEPIRWLLAAAGVEFEKFLKTRDDLARLRSDGSLMFQQVPMVEIDGMKLVQTKA
 ILNYIASKYNLYGKDMKERAIIDMYTEGVADLEIMILYYPHPPEEKEASLAKIKEQTRNRYFPAFEKVL
 KSHGQDYL VGNRLSRADIALVELLYHVEELDPGVVDNFP LLKALRSRVSNLPTVKKFLQPGSQRPFD
 DAKCVESAKKIFS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001077353

ORF Size: 1467 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001077353.2](#), [NP_001070821.1](#)

RefSeq Size: 1450 bp

RefSeq ORF: 666 bp


Locus ID: 14859

UniProt ID: [P30115](#)

Cytogenetics: 1 6.5 cM

Gene Summary: Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. This GST has a high catalytic activity for aflatoxin B1-8,9 epoxide. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG202521