

## Product datasheet for **MG202560**

### **Gsta2 (NM\_008182) Mouse Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Gsta2 (NM_008182) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Gsta2
Synonyms:	Gst2-2; Gstc-2; Gstc2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG202560 representing NM_008182 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCGGGAAGCCCGTGCTTCACCACTTCAATGCCCGGGCAGAATGGAGTGCATCAGGTGGCTCCTGG  
CTGCAGCAGGGGTGGAGTTTGAAGAGAAGTTTATACAGAGTCCGGAAGATTTGAAAAGCTAAAAAAGA  
TGGGAATTTGATGTTTGACCAAGTGCCCATGGTGGAGATTGATGGGATGAAGCTGGTGCAGACCAGAGCC  
ATTCTCAACTACATCGCCACCAATATGACCTCTATGGGAAGGACATGAAGGAGAGGCCCTGATTGACA  
TGTATTCAGAAGTATTTAGATCTGACTGAAATGATTGGCAATTGTTAATATGTCCCCAGACCAAAA  
AGAAGCCAAGACTGCCTTGGCAAAAGACAGGACCAAAAACCGTTACTTGCCTGCCTTTGAAAAGGTGTTG  
AAGGGCCATGGACAAGACTACCTTGTGGCAACAGGCTGACCAGGGTGGACGTCCACCTGCTGGAAGTTC  
TCCTCTATGTTGAAGAGTTTGTGCCAGCCTTCTGACTCCTTTCCTCTGCTGAAGGCCTTCAAGAGCAG  
AATCAGCAGCCTCCCAATGTGAAGAAGTTCCTACAGCCTGGCAGCCAGAGAAAGCCTCCCTTGGATGCA  
AAACAAATTGAAGAAGCAAGGAAGGTTTTCAAGTTT

**ACGCGT**ACGCGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG202560 representing NM\_008182  
 Red=Cloning site Green=Tags(s)

MAGKPVLHHFNARGRMECIRWLLAAAGVEFEFKFIQSPEDLEKLLKKGDNLMFDQVPMVEIDGMKLVQTRA  
 ILNYIATKYDLYGKDMKERALIDMYSEGILDLTEMIGQLLICPPDQKEAKTALAKDRTKNRYLPFAFEKVL  
 KGHGQDYLVGNRLTRVDVHLELLELLLYVEEFDASLLTPFPLLKAFKSRISLSPNVKFLQPGSQRKPPPLDA  
 KQIEEARKVFKF


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_008182

**ORF Size:** 666 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_008182.3</a> , <a href="#">NP_032208.2</a>
<b>RefSeq Size:</b>	855 bp
<b>RefSeq ORF:</b>	669 bp
<b>Locus ID:</b>	14858
<b>UniProt ID:</b>	<a href="#">P10648</a>
<b>Cytogenetics:</b>	9 43.65 cM
<b>Gene Summary:</b>	Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.[UniProtKB/Swiss-Prot Function]