

## Product datasheet for **MG202609**

### **Epdr1 (NM\_134065) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Epdr1 (NM_134065) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Epdr1
Synonyms:	AU040950; AV006210; Epdm2; Epdr2; MERP-1; MERP-2; MERP2; Ucc1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG202609 representing NM_134065 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCCGCGCGCTCCCCGCCGCTGGTCCAGGGGCTCGGGGACCTGGCTGCTGGGAAGCCTCTGGG  
TCTGGGTGCTGTGCGCCTGGGGATGGCGGGCTCCCTGGGAACCCACAGCCATGCCAGGCACCCAGCA  
GTGGGAGGGACGCCAGGTTCTGTACCAGCAGAGCAGCGGGCACAACAACCGCCCTGGTGTCTACGAT  
GGTCTCAACCAGCGCTGCGGGTCTGGACGAGAGGAAAGCGCTGATCCCCTGCAAGAGATTATTTGAAT  
ACATTTTACTCTATAAGGAGGGAGTGATGTTTCAGATTGAACAAGCCACCAACAGTGTGCAAAGATCCC  
CTTGGTGAATCCTGGGATCCTCTGGACATTCCCAGAATTCTACCTTTGAAGATCAGTACTCCATCGGA  
GGGCTCAGGAGCAGATCCTGGTCCAGGAGTGGTCTGACAGAAGAACAGCAAGATCCTATGAAACTTGG  
TCGGCGTTTATACAGCCAAGGATTGTTATCCGGTCCAGGAGACCTTCATCAGGAACACTGTTGGTTCAT  
GTCCACGCGGTTCTTTGATGTGCAGCTAGGCATTAAGGACCCCTCTGTGTTACCCACCAAGCAGCATGC  
CAGGCAGCGCAGCCAGAGAAGATGAGTGACGGCTGCTCCTTG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG202609 representing NM\_134065  
 Red=Cloning site Green=Tags(s)

MPARAPRRLVQGRGTWLLGSLVWVLCGLGMAGSLGTPQPCQAPQQWEGRQVL YQQSSGHNNRALVSYD  
 GLNQRVRLDERKALIPCKRLFYILLYKEGVMFQIEQATKQCAKIPLVESWDPLDIPQNSTFEDQYSIG  
 GPQEQIILVQEWSDRRTARSYETWIGVYTAKDCYPVQETFIRNYTVVMSTRFFDVQLGIKDPVFTFPSTC  
 QAAQPEKMSDGC SL

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_134065

**ORF Size:** 672 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_134065.4](#)

**RefSeq Size:** 2364 bp

**RefSeq ORF:** 675 bp


**Locus ID:** 105298

**UniProt ID:** [Q99M71](#)

**Cytogenetics:** 13 A2

**Gene Summary:** Binds anionic lipids and gangliosides at acidic pH.[UniProtKB/Swiss-Prot Function]

### Product images:



Circular map for MG202609