

Product datasheet for **MG202665**

Spcs2 (NM_025668) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spcs2 (NM_025668) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Spcs2
Synonyms:	5730406I15Rik; AA408713; mKIAA0102
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG202665 representing NM_025668 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGCGTCGGCGTCGCAGGGTGGGAGAAGCGCGGTGGCGGAGGCAGTAGTGGGGCTGGCGCGGCC
CCAGCTGCGGGACAAGCAGCAGCCGAGTGGTTTGTGGATAAGTGGAAGATTGATGATAAGCCTGTAAA
AATTGATAAGTGGGATGGATCAGCTGTGAAGAACTCTCTGGATGACTCTGCCAAAAGGTAAGTCTGGAA
AAGTACAAGTATGTGGAGAACTTTGGTCTCATTGACGGTCGCCTACCATCTGTACCATCTCCTGCTTTT
TTGCCATCGTGGCTTTGATTTGGGATTACATGCACCCCTTTCCTGAGCCCAAGCCAGTTTGGCTTTGTG
CGTCATATCCTATTTTGTGATGATGGGGATTCTGACCATTATACCTCATATAAAGAGAAGAGCATCTTC
CTTGTGGCCACAGAAAAGACCCACAGGAATGGATCCTGATGATATTTGGCAGCTGTCTCCAGCCTCA
AAAGGTTTGTGACAAATACACCCTGAAGTTGACCTTCATCAGTGGGAGGACAAAGCAGCAGCGGGAAGC
CGAGTTCACCAAGTCCATTGCTAAGTTCTTTGACCACAGCGGGACATTGGTCATGGATGCATATGAACCT
GAAATATCCAGGCTCCATGACAGTCTCGCCACAGAAAAGAAAAATAAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG202665 representing NM_025668
 Red=Cloning site Green=Tags(s)

MAASASQGGRRSGGGGSSGAGGGPSCGTSSSRSGLLDKWKIDDKPKIDKWDGSAVKNSLDDSAKKVLLE
 KYKYVENFGLIDGRLTICTISCFFAIVALIWDYMHPFPEPKPVLALCVISYFVMMGILTIYTSYKEKSIF
 LVAHRKDPTGMDPDDIWQLSSSLKRFDDKYTLKLTFFISGRTKQQREAEFTKSIAKFFDHSGLVMDAYEP
 EISRLHDSLATERKIK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_025668

ORF Size: 678 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_025668.2](#), [NP_079944.1](#)

RefSeq Size: 2816 bp

RefSeq ORF: 681 bp


Locus ID: 66624

UniProt ID: [Q9CYN2](#)

Cytogenetics: 7 E2

Gene Summary: Component of the microsomal signal peptidase complex which removes signal peptides from nascent proteins as they are translocated into the lumen of the endoplasmic reticulum. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG202665