

## Product datasheet for **MG202734**

### Rpe (NM\_025683) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Tag:** TurboGFP  
**Symbol:** Rpe  
**Synonyms:** 2810429B02Rik; 5730518J08Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)

**ORF Nucleotide Sequence:** >MG202734 representing NM\_025683  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGCGTCGGGCTGCAAGATCGGCCGTCCATCCTCAACAGCGACCTGGCCAACCTGGGGGCCGAGTGCC  
TGCGCATGTGGACTCGGGGGCCGACTACCTGCACCTGGATGTAATGGACGGGCATTTTGTCCCAATAT  
CACTTTTGGTCACCCTGTGGTTGAAAGCCTCCGAAAGCAGCTAGGCCAGGACCCTTTCTTTGATATGCAC  
ATGATGGTGTCTAGGCCAGAACAGTGGTAAAACCAATGGCTGTGGCAGGTGCCAATCAATATACCTTTC  
ATCTTGAGGCAACTGAGAACCAGGGCTTTGATTAAGACATCCGGGAGAATGGGATGAAGGTTGGCCT  
TGCCATCAAACCAGGAACACAGTTGAATTTGGCCCATGGGCTAATCAAATAGATATGGCCTTGGA  
ATGACCGTAGAACCTGGTGGGAGGGCAGAAATTCATGGAAGACATGATGCCAAGGTTCACTGGTTGA  
GGACCCAGTTCCAACCTTTGGATATAGAGGTCGATGGCGGAGTAGGTCCTGACACTGTTCAGAAGTGTG  
GGAGGCAGGAGCTAACATGATTGTGTCTGGCAGTGCCATTATGAGGAGTGATGACCCAGAGCTGTTATC  
AACTTACTAAGAAACGTTTGCTCAGAAGCGGCTCAGAAACGTTCCCTTGATCGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG202734 representing NM\_025683  
Red=Cloning site Green=Tags(s)


MASGCKIGPSILNSDLANLGAECLRMLDSGADYLHLDVMDGHFVFNITFGHPVVESLRKQLGQDPFFDMH  
MMVSRPEQWVKPMAVAGANQYTFHLEATENPGALIKDIRENGMKVGLAIKPGTTVEYLAPWANQIDMALV  
MTVEPGFGGQKFMEDMMPKVHWRQTQFPTLDIEVDGGVGPDTVQKCAEAGANMIVSGSAIMRSDDPRAVI  
NLLRNVCSEAAQKRLDR

TRTRPLE - GFP Tag - V



**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_025683

**ORF Size:** 684 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_025683.3](#)

RefSeq Size:	2631 bp
RefSeq ORF:	687 bp
Locus ID:	66646
UniProt ID:	<u>Q8VEE0</u>
Cytogenetics:	1 C3
Gene Summary:	Catalyzes the reversible epimerization of D-ribulose 5-phosphate to D-xylulose 5-phosphate. [UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MG202734