

## Product datasheet for **MG202754**

### Clec1b (NM\_019985) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Clec1b (NM_019985) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Clec1b
Synonyms:	1810061113Rik; Clec-2; Clec2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG202754 representing NM_019985 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCAGGATGAAGATGGGTATATCACTTTAAACATCAAGCCCCGAAACAAGCTCTCAGCTCAGCGGAAC  
CTGCCTCTTCTTGGTGGCGTGTGATGGCTTTAGTTCTGCTGATCTCATCCATGGGCTGGTTGTTGGACT  
CGTGGCTCTGGGGATCATGTCGGTCACACAGCAAAAGTATCTACTGGCGGAGAAGGAAAATCTCTCAGCG  
ACTCTGCAACAATTGGCCAAGAAATTCGCCAAGAGTTGATTAGACAATCAGAAATTAAGACAAAGAGCA  
CTTTTGAGCACAAGTGCAGCCCCGCGCCACGAAGTGGAGATACCATGGAGATAGTTGCTACGGGTTCTT  
CAGGCGTAACCTAACATGGGAAGAGAGCAAGCAGTATTGCACTGAGCAGAATGCAACACTTGTGAAGACT  
GCCAGCCAGAGCACCCCTGGACTACATTGCAGAAAGGATTACTTCAGTCCGTTGGATTGGATTATCAGGCC  
AGAACTCTAAGAAAGACTGGATGTGGGAGGATAGCTCAGTTCTTCGCAAGAACGGGATTAATCTTTCTGG  
GAATACAGAAGAAAACATGAATTGTGCTTATCTTCATAATGGAAAAATCCATCCAGCTTCTGTAAAGAG  
AGACATTACTTAATATGTGAGAGAAATGCTGGCATGACAAGAGTGGACCAACTGCTT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG202754 representing NM\_019985  
Red=Cloning site Green=Tags(s)

MQDEDEGYITLNIKPRKQALSSAEPASSWWRVMAVLLISSMGLVVGLVALGIMSVTQQKYLLAEKENLSA  
 TLQQLAKKFCQELIRQSEIKTKSTFEHKCSPCATKWRYHGDSCYGFRRNLTWEEKYQYCTEQNATLVKT  
 ASQSTLDYIAERITSVRWIGLSRQNSKKDWMWEDSSVLKNGINLSGNTENMNCAYLHNGKIHPASCKE  
 RHYLICERNAGMTRVDQLL

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_019985

**ORF Size:** 687 bp


**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_019985.3</a></u>
<b>RefSeq Size:</b>	1134 bp
<b>RefSeq ORF:</b>	690 bp
<b>Locus ID:</b>	56760
<b>UniProt ID:</b>	<u><a href="#">Q9JL99</a></u>
<b>Cytogenetics:</b>	6 F3
<b>Gene Summary:</b>	C-type lectin-like receptor that functions as a platelet receptor for the lymphatic endothelial marker, PDPN. After ligand activation, signals via sequential activation of SRC and SYK tyrosine kinases leading to activation of PLCG2.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MG202754