

Product datasheet for **MG202764**

Rnf141 (NM_025999) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rnf141 (NM_025999) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Rnf141
Synonyms:	2610110L04Rik; AA792898; AU022812; ZFP36; ZNF230
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG202764 representing NM_025999 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGACAGCAAATTTCCGATCAGACGCAGTTGGTTATTAACAAGTTACCAGAAAAGGTAGCAAAGCAGC
TCACCTTGGTTCGAGAAAGTGGCTCCTTAACCTACGAAGAGTTTCTTGGGAGAGTAGCTGAACCTCAATGA
TGTAACGTCTAAGGTGGCTGCTGGCCAGGAAAAGCATCTTCTTTGAGGTACAACCTGGGTCTGACTCC
TCTGCCTTCTGGAAAGTAGTTGTACGGTGGTCTGTACCAAGATTAACAAAAGCAGTGGCATTGTGGAAG
CATCGCGGATCATGAATTTGTACCAAGTTTATTCAGCTCTATAAAGACATCACAAGCCAAGCTGCTGGAGT
GCTGGCTCAGAGCTCCACCTCTGAAGAGCCTGATGAGAACCCATCCTCGGTCACATCTTGCCAGGCTAGT
CTATGGATGGGCAGGGTGAAGCAGCTGACTGATGAGGAGGAGTGTGTATCTGCATGGACGGCGAGCTG
ACCTCATCCTGCCCTGTGCTCACAGCTTCTGTGAGAAGTGCATTGATAAATGGAGTGATCGACACAGGAA
CTGCCCCATTTGTCGCCTACAGATGACTGGAGCAAATGAATCATGGGTGGTGTGATGACACCCACTGAA
GATGACATGGCTAACTATATTCTTAACATGGCTGATGAGGCAGGCCAGCCCCACAGGCCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG202764 representing NM_025999
 Red=Cloning site Green=Tags(s)

MGQQISDQTLVINKLPEKVAKHVTLVRESGSLTYEEFLGRVAELNDVTAKVAAGQEKHLLFEVQPGSDS
 SAFWKVVVRVCTKINKSSGIVEASRIMNLYQFIQLYKDITSQAAGVLAQSSTSEEPDENPSSVTSCQAS
 LWMGRVKQLTDEEECCICMDGRADLILPCAHSFQKCIDKWSDRHRNCPICRLQMTGANESWVSDAPTE
 DDMANYILNMADEAGQPHRP

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_025999

ORF Size: 690 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_025999.3](#)

RefSeq Size: 4019 bp

RefSeq ORF: 693 bp


Locus ID: 67150

UniProt ID: [Q99MB7](#)

Cytogenetics: 7 E3

Gene Summary: May be involved in spermatogenesis.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG202764