

Product datasheet for **MG202813**

Tk1 (BC012215) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Tk1 (BC012215) Mouse Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | Tk1 |
| Synonyms: | Tk1a, Tk1b |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >MG202813 representing BC012215 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCTACATCAATCTGCCACCGTGCTGCCAGCTCCCCAGCAAGACTCGGGGCAGATTCAGGTGA
TTCTCGGGCCCATGTTCTCAGGGAAAAGCACAGAGCTGATGAGAAGAGTCCGGCGCTTCCAGATCGCCCA
GTACAAGTGCCTGGTCATCAAGTATGCCAAGACACGCGCTATAGCAACAGCTTCTCCACACATGATCGG
AACACCATGGACGCATTGCCAGCCTGCATGCTCCGCGATGTGACCCAGGAGTCTTGGGTGTGGCCGTCA
TTGGCATCGATGAGGGGCAGTTTTTTTCTGACATTGTGGATTTCTGTGAAATGATGGCCAACGAGGGCAA
GACAGTAATTGTGGCAGCGCTGGATGGACCTTCCAGAGGAAGGCTTTCGGCAGCATCTTGAACCTGGTG
CCCCTGGCGGAGAGTGTGGTGAAGCTCACCGCTGTGTGCATGGAGTGCTTCCGAGAAGCTGCCTACACGA
AGAGGCTGGGCCTGGAGAAAGAGGTGGAGGTGATTGGCGGAGCCGACAAGTATCACTCCGTGTGCCCGCT
GTGCTACTTTAAGAAGCTTTCAGCCCAGACTGCTGGCTCAGACAACAAGAAGTGTCTGGTGTGGGGCAG
CCGGGAGAGGCCTTGGTTGTGAGGAAGCTTTTGCCTCTCAGCAAGTCTACAATACAACCTGCCAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG202813 representing BC012215
 Red=Cloning site Green=Tags(s)

MSYINLPTVLPSSPSKTRGQIQVILGPMFSGKSTELMRRVRRFQIAQYKCLVIKYAKDTRYSNSFSTHDR
 NTMDALPACMLRDVTQESLGVAVIGIDEGQFFPDI VDFCEMMANEGKTVIVAALDGTFRKAFGSILNLV
 PLAESVVKLTAVCMECFREAAAYTKRLGLEKEVEVIGGADKYHSVCRLCYFKKSSAQTAGSDNKNCLVLGQ
 PGEALVVRKLFASQQVLQYNSAN


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: BC012215

ORF Size: 701 bp

| | |
|-------------------------------|---|
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | BC012215 , AAH12215 |
| RefSeq Size: | 1175 bp |
| RefSeq ORF: | 701 bp |
| Locus ID: | 21877 |
| Cytogenetics: | 11 82.96 cM |