

Product datasheet for **MG202867**

Shisa5 (NM_025858) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Shisa5 (NM_025858) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Shisa5
Synonyms: 2310008D10Rik; 6430628I05Rik; mShisa5; Scotin
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG202867 representing NM_025858
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTGCGCCGGCGCCCTCTCTGTGGACCCTATTGCTGCTGCTGTTGCTGCTGCCGCCCTCCGGGTG
 CCCATGGTGAGCTGTGCAGGCCCTTTGGTGAAGACAATTCGATCCCAGTGTCTGTCTGATTTCTGTTG
 TGGTTCCTGTTCCAACCAATACTGCTGCTCGGACGTGCTGAGGAAAATCCAGTGAATGAGGAAATGTGT
 CCTGAGCCAGAGTCCAGCAGATTTCCACCCCGCGGAGGAGACACCCGAACATCTGGGTTACAGCGCTGA
 AATTTTCGATCCAGTTTTGACAGTGACCCTATGTCAGGGTTCGGAGCGACCGTCGCCATTGGCGTGACCAT
 CTTTGTGGTGTATTGCCACTATCATCATCTGCTTACCTGCTCCTGCTGCTGTGTATAAGATGTGC
 TGCCCCAACGCCCTGTCGTGACCAACACCACAACACTACTACCGTGGTTCATGCCCTTACCTCAGCCTC
 AACCTCAACCTGTGGCCCCAGCTATCCTGGACCAACATACCAGGGCTACCATCCCATGCCCCCCAGCC
 AGGAATGCCAGCAGCACCTACCCAACGAGTACCCACCACCTACCTGGCCCAGCCCACAGGGCCGCCA
 CCTACCATGAGTCCTTGGCTGGAGCCAGCCAGCCTCCATACAACCCGACCTACATGGATTCCCTAAAGA
 CAATCCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG202867 representing NM_025858
Red=Cloning site Green=Tags(s)

```


MAAPAPSLWLLLLLLLLPPPPGAHGLCRPFGEDNSIPVFCPDFCCGSCSNQYCCSDVLRKIQWNEEMC
PEPESSRFSTPAEETPEHLGSAKFRSSFDSDPMSGFGATVAIGVTIFVVFATIIICFTCCSCCLYKMC
CPQRPVVNTTTTTVVHAPYPQPQPVPVAPSYPGPTYQGYHPMPQPQGMPAAPYPTQYPPPYLAQPTGPP
PYHESLAGASQPPYNPTYMDSLKTIP
    
```

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_025858

ORF Size: 705 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_025858.2](#), [NP_080134.1](#)

RefSeq Size: 1792 bp

RefSeq ORF: 708 bp

Locus ID: 66940

UniProt ID: [Q9D7I0](#)

Cytogenetics: 9 F2

Gene Summary: Can induce apoptosis in a caspase-dependent manner and plays a role in p53/TP53-dependent apoptosis.[UniProtKB/Swiss-Prot Function]