

Product datasheet for **MG202900**

Dhrs7c (NM_001013013) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dhrs7c (NM_001013013) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Dhrs7c
Synonyms:	1110001P11Rik; AI120487; SRP-35
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG202900 representing NM_001013013 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGGCTCATGGCTGTCCTGATGCTACCCCTGCTGCTGCTGGGAATCAGCGGCCTCCTTTCATTTACC
AGGAGGCATCCAGGCTGTGGTCGAAGTCTGCCGTGCAGAACAAAGTGGTGGTCATCACAGATGCCATCTC
AGGACTGGAAAGGAGTGTGCTCGGGTGTCCATGCAGGTGGGCAAGGCTGGTGTGTGAAAGAAC
TGGGAGGGACTGGAGAGCCTCTATGCCACCTTGACCAGTGTGGCTGACCCAGCAAGACATTCAACCCCA
AGCTGGTCTCCTGGATCTCTCAGACATTAGCTGTGTTCAAGATGTGGCCAAAGAGGTCTGGACTGCTA
CGCTGTGTGGACATCCTCATCAACAATGCCAGGTGAAAGTGAAGGGGCTGCCACAAGATTTCCCTG
GAGCTTGACAAAAAGATCATGGATGCCAACTACTTCGGACCCATCACTTTAACCAAAGTTCTGTTCCCA
ACATGATCTCCAGGAGAACAGGCCAGATTGTGTTAGTGAACAACATCCAAGCGAAGTTTGGAAATCCCGTT
CCGCACAGCTTATGCAGCTTCTAAGCATGCCGTGCATGGGCTTCTTTGACTGCCTCCGAGCCGAGGTTGAG
GAATACGATGTTGTGGTCAGCACCGTGAGCCCACTTTTCCGCTCCTACCGTGCCTCCCTGAGCAAA
GAACTGGGAGACATCCATTTGTAATTCTTCTGCAGGAAGCTAGCCTATGGCGTGCACCCGGTGGAGGT
GGCTGAGGAAGTATGCGCACAGTACGGAGGAAGAAGCAAGAGGTGTTTCATGGCCAACCCGGTCCCAAG
GCTGCCGTGTTTCATCCGACCTTCTCCCTGAGTCTTCTTCGCTGTGGTGGCCTGTGGGTGAAGGAGA
AGCTCAATGTCCCAGAAGAGGGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG202900 representing NM_001013013
 Red=Cloning site Green=Tags(s)

MGLMAVLMPLLLLGISGLLFIYQEASRLWSKSAVQNKVVVITDAISGLGKECARVFHAGGARLVLCGKN
 WEGLESYATLTSVADPSKFTFPKLVLLDLSDISCVQDVAKEVLDCYGCVDILINNASVKVKGPAHKISL
 ELDDKIMDANYFGPITLTKVLLPNMISRRTGQIVLVNNIQAKFGIPFRTAYAASKHAVMGFFDCLRAEVE
 EYDVVVSTVSPTFIRSYRASPEQRNWETSICKFFCRKLAYGVHPVEVAEEVMRTVRRKKQEVFMANPVPK
 AAVFIRTFPEFFFAVACGVKEKLVNPEEG

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001013013

ORF Size: 933 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001013013.2](#), [NP_001013031.2](#)

RefSeq Size: 1129 bp

RefSeq ORF: 936 bp


Locus ID: 68460

UniProt ID: [Q8CHS7](#)

Cytogenetics: 11 B3

Gene Summary: Putative oxidoreductase.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG202900