

Product datasheet for **MG203032**

Schip1 (BC060529) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Schip1 (BC060529) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Schip1
Synonyms:	Nf2ip, SCHIP-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203032 representing BC060529 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTACATCAGGAGAACTGCTCGTACCAGGCACAGAAGAATGAGAGAGAGTCTATCAGACAGAAGTTGG
CACTCGGAAGCTTCTTTGACGATGGCCAGGAATCTATACCAGCTGCAGCAAAGTGGAAGCCAAGCCT
TTCTGCAAGACTACAGAGCGGGATGAACCTCCAGATATGCTTTGTCAATGACAGCGGCAGTGACAAGGAC
AGCGATGCAGATGACAGTAAGACGGAACACAGCTTGGACACGCCCTTGTCCTCCATGAGCAAGCAGAGTT
CTTCTATTCGGATAGAGACACAACCTGAGGAGGAGTCTGAATCCCTGGATGACATGGACTTCTCACAAG
GCAAAAGAAGCTACAAGCTGAAGCCAAAATGGCTCTGGCCATGGCCAAACCAATGGCCAAAATGCAAGTA
GAAGTGGAAAAACAGAACAGGAAAAAGTCTCCCGTCGCTGATCTTCTCCACACATGCCTCACATAAGCG
AATGTTTGATGAAAAGAAGCTTAAAGCCACCCGACCTGAGAGACATGACTATCGGGCAGCTACAAGTGAT
CGTCAATGACCTCCACTCCAGATTGAAAGTTTGAATGAAGAGTTGGTCCAGCTGCTCCTTATTCGAGAT
GAGCTGCACACAGAACAAGATGCCATGCTGGTGGACATTGAAGACTTGACTAGACACGCTGAGAGTCAGC
AGAAGCACATGGCTGAGAAAATGCCCGCAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203032 representing BC060529
 Red=Cloning site Green=Tags(s)

MVHQENCYSYAQKNERESIRQKLALGSSFFDDGPGIYTSCKSGKPSLSARLQSGMNLQICFVNDSGSDKD
 SDADDSKTETSLDTPSPMSKQSSYSRDTTEEESESLDDMDFLTRQKKLQAEAKMALAMAKPMAKMV
 EVEKQNRKSPVADLLPHMPHISECLMKRSLKPTDLRDMTIGQLQVIVNDLHSQIESLNEELVQLLLIRD
 ELHTEQDAMLVDIEDLTRHAESQQKHMAEKMPAK

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: BC060529

ORF Size: 734 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC060529 , AAH60529
RefSeq Size:	1569 bp
RefSeq ORF:	734 bp
Locus ID:	30953
Cytogenetics:	3 E1