

Product datasheet for **MG203038**

Ctdnep1 (NM_026017) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ctdnep1 (NM_026017) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ctdnep1
Synonyms:	2610507E10Rik; Dullard
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203038 representing NM_026017 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATCGGACGCAGTGTCTGCTGGGGCTGCGCACGTTTCGTGGCCTTCGCCGCAAGCTCTGGAGCTTCT
TTATTTACCTTTTTCGGAGGCAGATCCGCACGGTAATTCAGTATCAGACTGTTTCGATATGATATCCTGCC
CTTATCTCCTTTGTCCCAGGATCGCCTAGCCCAGGTGAAGAGGAAGATCCTGGTCTGGATCTGGACGAA
ACCCTGATTCACCTCTACCACGATGGGGTCTGAGGCCTACAGTGAGACCTGGGACACCTCCCGACTTCA
TCCTCAAGGTGGTAATAGACAAACACCCAGTCCGGTTTTTTGTACATAAGAGGCCCATGTGGATTTCTT
CTTAGAAGTGGTAAGCCAGTGGTATGAGCTTGTGGTGTTCACAGCAAGCATGGAAATTTATGGCTCTGCT
GTGGCAGATAAACTGGACAACAGCAGAAGCATTCTTAAGAGGAGATACTACAGACAGCACTGCACCTTTGG
AGTTGGGCAGCTACATCAAAGACCTCTCCGTGGTCCACAGCGACCTGTCCAGCATCGTGATCCTGGACAA
CTCCCCGGGGCTTACAGGAGCCACCCAGACAATGCCATCCCCATCAAATCCTGGTTCAGTGACCCAGT
GACACAGCCCTTCTCAACCTTCTCCAATGCTGGATGCCCTCAGGTTCACTGCTGATGTCGGATCGGTG
TGAGCCGAAACCTTACCAACATAGGCTCTGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203038 representing NM_026017
 Red=Cloning site Green=Tags(s)

MMRTQCLLGLRTFVAFAAKLWSFFIYLLRRQIRTVIQYQTVRYDILPLSPLSRNRLAQVKRKILVLDLDE
 TLIHSHHDGVL RPTV RPTGPPDFILKVVIDKHPVRFVHKRPHVDFLEVVVSQWYELVVFTASMEIYGSA
 VADKLDNSRSILKRYYRQHCTLELGSYIKDLSVVHSDLSSIVILDNSPGAYRSHPDNAIPIKSWFSDPS
 DTALLNLLPMLDALRFTADVRSVLSRNLHQHRLW

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_026017

ORF Size: 732 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_026017.2 , NP_080293.1
RefSeq Size:	1615 bp
RefSeq ORF:	735 bp
Locus ID:	67181
UniProt ID:	Q3TP92
Cytogenetics:	11 B3
Gene Summary:	Serine/threonine protein phosphatase forming with CNEP1R1 an active phosphatase complex that dephosphorylates and may activate LPIN1 and LPIN2. LPIN1 and LPIN2 are phosphatidate phosphatases that catalyze the conversion of phosphatidic acid to diacylglycerol and control the metabolism of fatty acids at different levels. May indirectly modulate the lipid composition of nuclear and/or endoplasmic reticulum membranes and be required for proper nuclear membrane morphology and/or dynamics. May also indirectly regulate the production of lipid droplets and triacylglycerol. May antagonize BMP signaling (By similarity).[UniProtKB/Swiss-Prot Function]