

## Product datasheet for **MG203046**

### **Pdc (NM\_024458) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Pdc (NM\_024458) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Pdc  
**Synonyms:** Rpr-1; Rpr1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG203046 representing NM\_024458  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAAGAAGCCGCCAGCCAAAGCTTAGAAGAAGATTTTGAAGGACAGGCCACGCACACAGGACCCAAAG  
GAGTAATAAATGATTGGAGAAAGTTTAAATTAGAAAGTGAAGATGGTGATTCAATCCACCCAGCAAGAA  
GGAAATCCTCAGACAAATGTCCTCTCCTCAGAGCAGAGATGACTCAAAGAAAGAATGAGCAGAAAGATG  
AGCATTCAAGAATATGAACTCATTCAAGACAAGAAGATGAGAGTTGCCTTCGCAAATACCGGAGGC  
AGTGCATGCAGGATATGCATCAGAAGTTGAGTTTGGGCCTAGGTATGGGTTTGTGTATGAGCTGGAAC  
AGGGGAGCAATTCCTGGAAACATCGAGAAGGAGCAGAAGGTACCACAATTGTGGTTAACATTTACGAG  
GATGGTGTGAGGGCTGCGATGCACTCAACAGCAGCTTAGCGTGCCTTGCAGTAGAGTACCAATGGTCA  
AGTTCTGTAAAATAAAAGCTTCTAATACCGGAGCGGGGACCGCTTTTCCACAGACGTA CTCCCTACGTT  
GCTGGTGTACAAAGGTGGGAACTCATAAGCAATTTTATTAGTGTGGCTGAACAATTTGCTGAAGAATTT  
TTTGCTGTGGATGTGGAGTCTTCTAAATGAATATGGATTACTACCAGAAAGAGAGATACATGACCTAG  
AGCAGACCAACATGGAAGACGAGGACATTGAA

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG203046 representing NM\_024458  
Red=Cloning site Green=Tags(s)

MEEAASQSLLEEDFEGQATHTGPKGVINDWRKFKLESEDGDSIPPSKKEILRQMSSPQSRDDSKERMSRKM  
 SIQEYELIHQDKEDESLRKYRRQCMQDMHQKLSFGPRYGFVYELETGEQFLETIEKEQKVTIIVVNIYE  
 DGVRGCDALNSSLACLAVEYPMVKFCKIKASNTGAGDRFSTDVLPDLLVYKGGELISNFSVAEQFAEEF  
 FAVDVESFLNEYGLLPEREIHDLLEQTNMEDEDIE

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_024458

**ORF Size:** 1248 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_024458.2</a> , <a href="#">NP_077778.1</a>
<b>RefSeq Size:</b>	1263 bp
<b>RefSeq ORF:</b>	735 bp
<b>Locus ID:</b>	20028
<b>UniProt ID:</b>	<a href="#">Q9QW08</a>
<b>Cytogenetics:</b>	1 63.84 cM
<b>Gene Summary:</b>	Inhibits the transcriptional activation activity of the cone-rod homeobox CRX (By similarity). May participate in the regulation of visual phototransduction or in the integration of photoreceptor metabolism.[UniProtKB/Swiss-Prot Function]