

## Product datasheet for **MG203063**

### Tspan6 (NM\_019656) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Tspan6 (NM_019656) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Tspan6
Synonyms:	6720473L21Rik; AI316786; T245; Tm4sf; Tm4sf6; TSPAN-6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203063 representing NM_019656 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGTCCCATCTCGGAGGCTGCAGACTAAGCCGGTCATTACTTGTGGAAAAGCGTTCTCTTGATCT  
ACACTTTCATCTTTGGATCACTGGTGTATCCTTCTTGCCGTTGGCATTGGGGCAAGGTGAGCCTGGA  
AAATTATTTTCCCTCTAAATGAGAAGGCCACCAATGTCCCTTTGTGCTCATTGGCACAGGCACTGTC  
ATTATTCTCTGGGCACCTTTGGTTGTTTCGCTACCTGTCGAGCCTCTGCCTGGATGCTGAAACTGTACG  
CGATGTTTCTGACTCATTTTTTGGTGAAGTACTGTCGCGCCATTGTTGGATTGTTTTAGACATGA  
GATTAAGAAGCAGCTTTAAAAGTAACTATGAAAACGCTCTGAAGGAGTACAACCTCCACAGGAGACTATAGA  
AGTGAAGCTGTAGACAAGATCCAAAGTACGTTGCATTGTTGCGGTGTACCAATTACGGAGATTGAAAAG  
GTACGAACTATTACTCAGAAACAGGATTTCCCAAGAGCTGCTGTAAGTGGAGGGCTGTTATCCACAGAG  
AGATGCAGATAAAGTCAACGAAGAGGGTTGTTTCATAAAGGTAATGACAACCTATAGAGTCAGAAATGGGA  
GTCGTTGCTGGAATTTCTTTGGAGTTGCTTTCAGCTAATTGGAATCTTTCTTGCTACTGCCTGT  
CTCGTGCCATAACGAATAACAGTATGAGATAGTA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203063 representing NM\_019656  
 Red=Cloning site Green=Tags(s)

MASPSRRLQTKPVITCLKSVLLIYTFIFWITGVILLAVGIWGVSLNENYFSLNEKATNVPFVLIGTGTV  
 IILLGTFGCFATCRASAWMLKLYAMFLTLIFLVELVAAIVGVFRHEIKNSFKSNYENALKEYNSTGDYR  
 SEAVDKIQSTLHCCGVTNYGDWKGNTNYYSETGFPKSCCKLEGCYPQRDADKVNNEEGCFIKVMTTIESEMG  
 VVAGISFGVACFQLIGIFLAYCLSRAITNNQYEIV

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM\_019656

ORF Size: 735 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_019656.2</a>
<b>RefSeq Size:</b>	1763 bp
<b>RefSeq ORF:</b>	738 bp
<b>Locus ID:</b>	56496
<b>Cytogenetics:</b>	X E3