

Product datasheet for **MG203160**

Uevld (BC054796) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Uevld (BC054796) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Uevld
Synonyms:	8430408E05Rik; Attp
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203160 representing BC054796 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTTCGACTGTGAGGGCGTGAGACGGCTGCTAGGCAAGTACAAGTTCAGGGATCTGACTGTGGAAG
AGCTGAATGTGAGCGTGTCTTTCCCGCACTTCAGATACTCCGTGGACACCTATGTTTTCAAAGATACATC
CCAAAAAGACCTGCTGAATTTACCGGTACAATTCCTGTGATGTATCAGGGTAAGACATACAACATACCC
ATTCGCTTCTGGATTTTGGATTTCCACCCCTTTGCTCCACCCATTGCTTCTTGAAGCCAACCTGCAACA
TGGAATCTCGGTTGAAAACATGTGGATGCCAAAGGCAGAATATACTTGCCTATCTCCAAAACCTGGAG
CCACCCCAAACTGCCATCGTTGGATTAATTAAGAAATGATTGCCAAGTTTCAAGAGGAACTTCCTCTC
TATTCATACCGTCGTCCAATGAAGCGCAACAGGTGGACTTGCTGGCCTATATCAGAAAAATCACGGAAG
GTGTCTCAGACATAAATTCAAGGGGCTGGACAAATCATGAGAATAAAATACTCAATAAAATCACTGTGGT
TGGAAGTGGAGATTTGGGAATCGCTGCACATTGGCAATTTAGCAAAAGGGCATTGCAGACAAGCTGCTC
CTCTTAGACCTCTCAGACGGGATGAGCCAAGGGACGATGGACCTTGATATCTTCAACCTGCCTAATGTAG
AGATCAGCAAAGGTGGTATCTTCACAGCCAACCTTTGGGCGGTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)


Protein Sequence: >MG203160 representing BC054796
 Red=Cloning site Green=Tags(s)

MEFDCEGVRLLGKYKFRDLTVEELNVSVSFPHFYSVDITYVFKDTSQKDLLNFTGTIPVMYQGKTYNIP
 IRFWILDSHPFAPPICFLKPTANMEISVGKHVDAGRIYLPYLQNSHPKSAIVGLIKEMIAKFQEELPL
 YSIPSSNEAQQVDLLAYITKITEGVSDINSGWTNHNENKILNKITVVGSGDLGIACTLAISAKGIADKLL
 LLDLSDGMSQGTMDLDIFNLPNVEISKGGDLHSQLFGRF

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC054796

ORF Size: 747 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC054796](#), [AAH54796](#)

RefSeq Size: 4460 bp


RefSeq ORF: 749 bp

Locus ID: 54122

Cytogenetics: 7 B3

Gene Summary: Possible negative regulator of polyubiquitination.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG203160