

Product datasheet for **MG203166**

Taf9b (BC066223) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Taf9b (BC066223) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Taf9b
Synonyms:	BC066223; Taf9l
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203166 representing BC066223 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGCCGGCCAAGATGGCGCCTATCAAGAACGCTCCGAGAGATGCCTTGGTGATGGCACAGATCCTGA
AGGATATGGGAACCACAGAGTATGAGCCAAGGTTATAAATCAAATGTTGGAATTTGCTTTCCGATATGT
GACTACAATTCTGGACGATGCAAAGATTTATTCGAGTCATGCTAAGAAACCCACTGTTGATGCCGATGAC
GTGAGACTGGCCATCCAGTGTGGGCTGACCAGTCCTTTACCTCTCTCCCAAGAGATTTTTACTAG
ATATTGCAAGGCAGAAGAATCAAACCCCTTTGCCGCTGATTAAGCCATATGCAGGACCCAGACTGCCACC
TGATAGATACTGCTTAACAGCTCAAACCTATAGGCTGAAGTCCTTAGTTAAAAAGGGACCTAACCAAGGA
AGACTAGTGCCTAGATTAAGTGGGTCAGTAGTAGACCACTACTCCTCTGTAGCACCCCAAGCAG
TATCTGGCCAAATAAGGCTGCAACTCCGGTGTGAGTACAAGCCAAAGATTTGCTGTGCAGATTCACC
TTCTCAGTCCACACCTGCCAAGCCAGCTCCTGCAGCACTGCAGTCCAAAATGTTCTGATTAATCCTTCA
ATGATTGGGCCAAAAATATTCTCATTACCACCAGCATGGTTTCCTCACAGAACACAGCCACCGACTCAA
ACCCACTGAAGAGGAAGCATGATGACGACGATGACAATGACACCATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203166 representing BC066223
Red=Cloning site Green=Tags(s)


MEPAKMAPIKNAPRDALVMAQILKDMGTTEYEPRVINQMLEFAFRYVTTILDDAKIYSSHAKKPTVDADD
 VRLAIQCRADQSFTSPPPRDFLLDIARQKNQTPPLPIKPYAGPRLPPDRYCLTAPNYRLKSLVKKGPNQG
 RLVPRLSAVSSRPTTPPVAPPQAVSGPNKAATPVSVTSQRFVAVQIPPSQSTPAKPAPAATAVQNVLINPS
 MIGPKNILITTSMVSSQNTATDSNPLKRRKHDDDDNDTM

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: BC066223

ORF Size: 749 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC066223 , AAH66223
RefSeq Size:	2393 bp
RefSeq ORF:	749 bp
Locus ID:	407786
Cytogenetics:	X D
Gene Summary:	Essential for cell viability. TAF9 and TAF9B are involved in transcriptional activation as well as repression of distinct but overlapping sets of genes. May have a role in gene regulation associated with apoptosis. TAFs are components of the transcription factor IID (TFIID) complex, the TBP-free TAFII complex (TFTC), the PCAF histone acetylase complex and the STAGA transcription coactivator-HAT complex. TFIID or TFTC are essential for the regulation of RNA polymerase II-mediated transcription (By similarity).[UniProtKB/Swiss-Prot Function]