

Product datasheet for **MG203184**

Cyb561 (NM_007805) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cyb561 (NM_007805) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Cyb561
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG203184 representing NM_007805
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGGATCGCC**

ATGGAGCACAGTTCTGCGTCTGTCCCTGCTGCACTGCCGTACTATGTGGCCTTCTCCAGCTGCTGGGCC
 TCACTGTAGTGGCTGTGACTGGTGCCTGGCTGGTCTGTACCGAGGTGGCATTGCCTGGGAAAGCTCCCT
 GCAGTTTAATGTGCACCCGCTCTGCATGGTCATAGGCATGATCTTCCTCAAGGAGATGCTCTGCTGGTG
 TACCGTGTCTTCAGAAGAGAAGCCAAACGCACGACCAAGATCTTGCATGGGCTGCTGCACGTCTTCGCCT
 TCATCATCGCCCTGGTGGGTTGGTGGCTGTGTTTCGACTACCACAAGAAGAAGGGCTATGCAGACCTGTA
 CAGCTTGACAGCTGGTGTGGGATCTTAGTCTTTGTTCTGTACTTTGTGAGTGGCTCGTGGGTTTCAGC
 TTCTTCTGTTCCCTGGAGCTTCATTCTCTGCGGAGCCGCTACCGCCCTCAGCACATTTCTTTGGTG
 CCACCATCTTCTCTCTGTGGGCACAGCCCTACTGGGCTGAAGGAGGCCCTGCTGTTAAATTGGG
 GTCCAAGTACAGCACGTTTGAACCTGAGGGGGTGTGGCAACGTGCTAGGCTTGCTGCTGGTCTGTTTT
 GGGGTGGTTGTCTCTACATCTTGGCTCAGGCTGACTGGAAGCGGCCCTCCAGGCTGAAGAGCAAGCCC
 TCTCCATGGACTTAAGACACTACGGAGGGAGACAGCCCCAGTCCCCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG203184 representing NM_007805
 Red=Cloning site Green=Tags(s)

MEHSSASVPAALPYVAFSQLGLTVVAVTGAWLGLYRGGIAWESSLQFNVHPLCMVIGMIFLQGDALLV
 YRVFRREAKRTTKILHGLLHVFAFIILVGLVAVFDYHKKGYADLYSLHSWCGILVFVLYFVQWLVGFS
 FFLFPGASFSLRSRYRPQHIFFGATIFLFSVGTALLGLKEALLFKLGSKYSTFEPEGVLANVLGLLLVCF
 GVVVLYILAQADWKRPQAEQALSMDFKLTTEGDSPPQ


TRTRPLE - GFP Tag - V



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_007805

ORF Size: 750 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_007805.5](#)

RefSeq Size: 2475 bp

RefSeq ORF: 753 bp

Locus ID: 13056
 UniProt ID: [Q60720](#)
 Cytogenetics: 11 68.81 cM
 Gene Summary: Secretory vesicle-specific electron transport protein.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG203184