

Product datasheet for **MG203186**

Kcnip4 (BC051130) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kcnip4 (BC051130) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Kcnip4
Synonyms:	KChIP4, Calp250, Calp
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203186 representing BC051130 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAACGTGAGAAGGGTGGAAAGCATCTCGGCTCAGCTGGAGGAGGCGAGCTCCACAGGCGGTTTCCTCT
ACGCGCAGAACAACCAAGCGCAGCATTAAAGAGCGGCTCATGAAGCTTTGCCCTGCTCAGCTGCCAA
AACGTCGCTCCTGCTATTCAAACAGCGTGGAAGATGAGCTGGAGATGGCTACTGTCAGGCATCGGCCT
GAAGCCCTGGAGCTGCTGGAGGCCAGAGCAAATTCACCAAGAAAGAGCTTCAGATTCTTACAGAGGAT
TTAAGAATGAATGCCCCAGTGGTGTGTTAATGAAGAACTTTCAAGGAGATTTACTCACAGTTCTTTCC
ACAGGGAGACTCCACCACATATGCACATTTTCTCTTCAATGCATTGACACGGACCAACAATGGAGCTGTG
AGCTTTGAGGATTTATCAAAGGTCTTTCCATTTTGCTTCGAGGGACAGTACAAGAAAACTCAACTGGG
CATTTAATTTGTATGACATAAACAAAGATGGCTACATCACTAAAGAAGAAATGCTGGACATAATGAAAGC
AATCTACGACATGATGGGAAATGCACATACCCGGTCTCAAGGAAGATGCTCCCCGACAGCATGTGGAG
ACGTTCTCCAGAAGATGGACAAAAATAAAGATGGTGTGTTACCATAGATGAGTTTCATTGAAAGTTGCC
AAAAAGATGAAAACATAATGCGCTCCATGCAGCTCTTTGAAAATGTGATC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203186 representing BC051130
Red=Cloning site Green=Tags(s)

MNVRVESISAQLEEASSTGGFLYAQNNTKRSIKERLMKLLPCSAAKTSSPAIQNSVEDELEMATVRHRP
 EALELLEAQSKFTKKELQILYRFGKNECPSGVVNEETFKEIYSQFFPQGDSTTYAHFLFNAFDTDHNGAV
 SFEDFIKGLSILLRGTVQEKLNWAFNLVDINKDGYITKEEMLDIMKAIYDMMGKCTYPVLKEDAPRQHVE
 TFFQKMDKNKDGVVTTIDEFIESCQKDENIMRSMQLFENVI

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: BC051130

ORF Size: 752 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC051130 , AAH51130
RefSeq Size:	2392 bp
RefSeq ORF:	752 bp
Locus ID:	80334
Cytogenetics:	5 B3
Gene Summary:	Regulatory subunit of Kv4/D (Shal)-type voltage-gated rapidly inactivating A-type potassium channels, such as KCND2/Kv4.2 and KCND3/Kv4.3 (PubMed:19109250). Modulates channel expression at the cell membrane, gating characteristics, inactivation kinetics and rate of recovery from inactivation in a calcium-dependent and isoform-specific manner. [UniProtKB/Swiss-Prot Function]