

## Product datasheet for **MG203213**

### 1700074P13Rik (NM\_028550) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	1700074P13Rik (NM_028550) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	1700074P13Rik
Synonyms:	AV047902; Tryx5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203213 representing NM_028550 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAAGACTTTTTTCATCTTCACTCTCTTAAGTCTCGCAGTAGCTAATACTCCCAAAGCCGTTCTGAAAG  
AATATAACAATTCAGAAGAATACTTACCAGAGAATTTAATGTTCTTATATGGTCTATCTGCAGTCCAG  
CCCAGAGCCTTGTGTTGGGACTCTCATTGACCCTCAGTGGGTGCTCACCGCTGCTCATTGTTCTTTACCA  
ACTAAAATCCGACTGGGAGTTTACCGACCCAACATCAAAAATGAGAAAGAACAGATATGTAATTACTCAT  
TCACTGTAGTCCATCCTAATTTTGATGCAAAATTAAGTAAAAATGACTTGATGCTGATAAACTCTCATA  
TCCTGCTACTATCAACATGTATGTGGAACTATAGCCATAGCTATGGAACCAATGGCATTAAATGAATCC  
TGCTTCATACCAACATGGACCTGGAATAACTACAAAACCTTAGTGATCCTGACATCTTGACATGGATAA  
ATGAGTATTCTCTTTCCCAAGTGAAGTGTGGGACTCTTCAACAGAAACAAGAAACAAGAAATAAA  
CATTATGTGCATAGGACATTCACCTAATGCTATGTCTGCAACTAAGGAAGTATCAGCTGCTCCAGCCATC  
TGCAGTGGGAGGGTGCATGGAATATTGCTCTGGGGAAAGCAAGTGTAGCCAATGGAAGCAAAGGCTTCT  
TCACTGAAATTCATCCCTATGCAAGATGGATCTTAAAAATGATACGATCACAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG203213 representing NM\_028550  
 Red=Cloning site Green=Tags(s)


MKTFFIFITLLSLAVANTPKAVLKEYNNSEEYLPENFNVPYMYVLQSSPEPCVGTLIDPQWVLTAAHCSLP  
 TKIRLGVYRPNIKNEKEQICNYSFTVVHPNFDAKLLKNDLMLIKLSYPATINMYVGTIAIAMEPMAFNES  
 CFIPTWTWNNYKNLSDPDILTWINEYSLSPSDCLDLHQKQKETRINIMCIGHSLNAMSATKEVSAAPAI  
 CSGRVHGILSWGKASVANGSKGFFTEIHPYARWILKMRSH

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_028550

**ORF Size:** 753 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_028550.4</a>
<b>RefSeq Size:</b>	1647 bp
<b>RefSeq ORF:</b>	756 bp
<b>Locus ID:</b>	73481
<b>Cytogenetics:</b>	6 B1