

Product datasheet for **MG203217**

2310026E23Rik (BC051512) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	2310026E23Rik (BC051512) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	2310026E23Rik
Synonyms:	2310026E23Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203217 representing BC051512 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACAGTCTGGAGGATGCTGTCGTGCCCGGGCTCTGTATGAGGAGCTGCTGCGAACTACCAGCAGC
AACAAAGAGGAGATGCGCCACCTCCAGCAGGAGCTGGAGCGGACTCGGAGGCAGCTCGTACAGCAGGCCAA
GAAGCTCAAGGAGTACGGGGCATTGGTGTCTGAAATGAAGGAGCTCCGAGACCTCAACAGGAGACTCCAA
GATGTGCTACTCCTGCGGCTTGGCAGCGGTCCCGCCATTGACTTGGAAAAAGTAAAGTCAGAATGTCTCG
AGCCCCGAGCCGGAGTTACGGAGCACTTTCAGTGAGGAAGCAAATACGTCGTCCTATTACCCCGCTCCTGC
GCCTGTGCATGGACAAGTATATCCTAGACAATGGCAAGGTCCATCTGGGAAGCGGGATTTGGGTTGATGAG
GAGAAATGGCACCAGCTACAAGTTACCCAAGGAGATTCTAAGTACACGAAGAAGTGGCAGTCATGATCT
GGGGAACAGATGTTCTGAAAAACAGAAGTGTACAGGAGTCGCCACAAAAAAGAAAGGATGCAATCCC
AAAGCCACCCCTCTCCTCACAAGCTCAGCATCGTCAGAGAGTGTATATGACAGAATAGCACAAAGAA
ACGGTGGATGAAACGAAATTGCACAGAGACTCTCAAAGTCAACAAGTACATCTGTGAGAAAATCATGG
ATATCAACAAATCTTGTA AAAATGAAGAACGAAGGAAGCAAAAATACAATTTGCAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203217 representing BC051512
 Red=Cloning site Green=Tags(s)

MDSLEDAVVPRALYEELLRNYQQQEEMRHLQQELERTRRQLVQQAQKLKEYGALVSEMKELRDLNRRLO
 DVLLLRGSGPAIDLEKVKSECLEPEPELRSTFSEEANTSSYYPAPAPVMDKYILDNGKVHLGSGIWDDE
 EKWHQLQVTQGD SKYTKNLAVMIWGTDLKNRSVTGVATKKKKDAIPKPPLSPHKLSIVRECLYDRIAQE
 TVDETEIAQRLSKVNKYICEKIMDINKSCKNEERREAKYNLQ

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: BC051512

ORF Size: 758 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC051512 , AAH51512
RefSeq Size:	930 bp
RefSeq ORF:	758 bp
Locus ID:	67621
Cytogenetics:	4 D1
Gene Summary:	Acts as a transcriptional repressor.[UniProtKB/Swiss-Prot Function]