

## Product datasheet for **MG203274**

### Lrrc18 (BC061072) Mouse Tagged ORF Clone

#### Product data:

|                           |                                                                            |
|---------------------------|----------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                        |
| Product Name:             | Lrrc18 (BC061072) Mouse Tagged ORF Clone                                   |
| Tag:                      | TurboGFP                                                                   |
| Symbol:                   | Lrrc18                                                                     |
| Synonyms:                 | 4930442L21Rik; MTLR1                                                       |
| Mammalian Cell Selection: | Neomycin                                                                   |
| Vector:                   | pCMV6-AC-GFP (PS100010)                                                    |
| E. coli Selection:        | Ampicillin (100 ug/mL)                                                     |
| ORF Nucleotide Sequence:  | >MG203274 representing BC061072<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCAAAGGAGGGAAAGGCCCAAGGGCAAAAAGATCACCTTAATGTGGCAAAGAACTGCATCAAAA  
TCACATTTGATGGGAGAAAACGCCTTGACCTGAGCAAGATGGGTATCACCACTTCCCAAGTGTATCTT  
GCGGCTCAGTGATATAGATGAGCTAGATCTTAGCCGGAATATGATCAGAAAAATCCCGACTCCATCGCC  
AAGTTCCAGAACCTGCGATGGCTCGACCTGCACAGTAATATTTGACAAGCTCCCTGAATCCATTGGCC  
AGATGACCTCCTTGCTCTTCTCAATGTGAGCAACAACAGGCTGACCACCAACGGGCTTCTGTGGAGCT  
GAACCAACTCAAGAACATCCGCACCGTGAATCTGGGCCTTAACCACCTGGACAGCGTGCCCAACCACTG  
GGGGCCCTGAAGGAGCTCCATGAGGTGGGGCTGCATGACAATTTGCTGACTACCATCCCCGCCAGCATCG  
CCAAGTCCCCAAGCTGAAGAAGCTCAACATAAAGAGGAACCCCTTCCCAATGCAGATGAATCGGAGAT  
GTTCGTAGACTCCATCAAAGGCTAGAAAACCTCTATCTGGTGAAGAGAAGGATATGTGTTTCATCCTGT  
CTGCAGAGATGCCAACAGGCCAGGGACAAGTTGAATAAAATCAAGAGCATGGCCCCCTCTGCACCGAGAA  
AGGCCCTCTTTTCCAATTTGGTTTACCCAACCTCAACAGCCAAGGATGCCAGGAAGAATGGAGG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)

**Protein Sequence:** >MG203274 representing BC061072  
 Red=Cloning site Green=Tags(s)

MAKGGKGPKGKKITLNVAKNCIKITFDGRKRLDSLKMGITTFPKCILRLSDIDELDL SRNMIRKIPDSIA  
 KFQNLRWLDLHSNYIDKLPESIGQMTSLLFLNVSNNRLTTNGLPVELNQLKNIRTVNLGLNHLDSVPTTL  
 GALKELHEVGLHDNLLTTIPASIAKLPLKLNKRNPFNADESEMFDVSIKRLNLYLVEEKDMCSCS  
 LQRCQQARDKLNKIKSMAPSAPRKALF SNLVSPNSTAKDAQEEWR

TRTRPLE – GFP Tag – V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** BC061072

**ORF Size:** 767 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC061072](#), [AAH61072](#)

**RefSeq Size:** 1007 bp


**RefSeq ORF:** 767 bp

**Locus ID:** 67580

**Cytogenetics:** 14 B

**Gene Summary:** May be involved in the regulation of spermatogenesis and sperm maturation.  
[UniProtKB/Swiss-Prot Function]

## Product images:



Circular map for MG203274