

Product datasheet for **MG203321**

Enoph1 (NM_026421) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Enoph1 (NM_026421) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Enoph1
Synonyms:	2310057D15Rik; BB183658; C81437
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203321 representing NM_026421 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTCGTGGTTCCGTGCCCGCCGAAGTCACCGTAATCCTGTTAGATATCGAAGGTACCACAACCCCGA
TTGCTTTCGTGAAGGATGTTTTATTTCCGTACATCAAAGAAAATGTCAAGGAGTACCTGCAGACCCACTG
GGAGGAGGAGGAGTGTGAGCAGGACATCAGCCTGCTGAGGAAACAGGCCGAGGAAGATGCACATCTGGAC
GGGGCTGTTCTGATACCCGTGGCTTCTGGCAGTGATCTGCAGCAGATGATCCAGGCTGTGGTGGACAACG
TGTACTGGCAGATGTCCCATGACCGAAAGACCACAGCGCTCAAACAGCTGCAAGGCCACATGTGGAAGGC
AGCTTTCACGGCTGGCCGCATGAAGGCAGAGTTCTTTGCAGATGTGGTTCCTGCAGTCAGAAGGTGGAGA
GAGGCTGGGATGAAGGTATACATCTATTCTCAGGGAGCGTGGAGGCACAGAAGCTGCTCTTCGGCCATT
CCACAGAGGGAGACATTCTTGAGCTCATCGACGGTCACCTTTCACCAAGATTGGACACAAGTGGACAG
TGAGAGTTACCGGAAGATTGCTGATAGCATCGGTTGTTCAACCAACAACATCCTGTTTCTCACGGATGTC
ACCGTAGAGGCCAGTGCTGCTGAGGAGCGGACGTGCACGTGGCTGTGGTGGTGGCAGCTGGCAACGCGG
GGCTCACTGACGACGAGAAGACATACTACAACCTCATCACGTCTTCAGCGAGCTCTACCTGCCGTCCAC
C

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203321 representing NM_026421
Red=Cloning site Green=Tags(s)

MVVVSVPAEVTVILLDIEGTTTPIAFVKDVLFPYIKENVKEYLQTHWEEECQQDISLLRKQAEEDAHLDF
 GAVLIPVAGSGLDQQMIQAVVDNVYQMSHDRKTTALKQLQGHMWKAAFTAGRMKAEFFADVVPVRRWR
 EAGMKVYIYSSGSVEAQKLLFGHSTEGDILELIDGHFDTKIGHKVDSESYRKIADSIGCSTNNILFLTDV
 TVEASAAEEADVHVAVVVRPGNAGLTDDEKTYNLTFSSELYLPST

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_026421

ORF Size: 771 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_026421.1 , NP_080697.1
RefSeq Size:	1809 bp
RefSeq ORF:	774 bp
Locus ID:	67870
UniProt ID:	Q8BGB7
Cytogenetics:	5 E4
Gene Summary:	Bifunctional enzyme that catalyzes the enolization of 2,3-diketo-5-methylthiopentyl-1-phosphate (DK-MTP-1-P) into the intermediate 2-hydroxy-3-keto-5-methylthiopentyl-1-phosphate (HK-MTPenyl-1-P), which is then dephosphorylated to form the acireductone 1,2-dihydroxy-3-keto-5-methylthiopentene (DHK-MTPene).[UniProtKB/Swiss-Prot Function]