

Product datasheet for **MG203397**

H2-DMa (NM_010386) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: H2-DMa (NM_010386) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: H2-DMa
Synonyms: H-2Ma; H2-Ma
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG203397 representing NM_010386
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGCATGAGCAGAAGTCAGGAGCTGTGCTGCTGAGGCTGCTGCGCCTTCTGTGGCTGCTGCCCCACT
CCTGGGCGGTGCTCGAAGCATCTACACCAGTGTCTGGGATGACCCACAGAACCACACATTCCGGCACAC
TCTATTCTGCCAGGACGGGATCCCAACATAGGGCTCTCGGAGACCTATGACGAGGACGAACCTTCTCC
TTCGACTTCTCCAGAACACCAGAGTGCCCGACTGCCTGACTTTCAGAGTGGGCTCAGGGACAGGGG
ATGCCTCTGCCATTGCATTTGACAAAAGCTTCTGCGAGATGTTGATGCGGGAAGTGAGCCGAAGCTTGA
AGGTCAAATCCCAGTGTCCAGAGGTTTGCTGTGGCAGAGGTATTCACACTGAAGCCCCTGGAGTTTGGC
AAGCCCAACACGTTGGTCTGTTTCATCAGCAACCTCTTTCCACCGACCCTGACCGTGAACGGCAGCTTC
ACTCGGCCCTGTGGAGGGAGCCAGCCCCACATCCATCTCAGCTGTGATGGGCTGACCTTCAGGCTTT
CTCTTATTTAAACTTCACACCGGAACCTTTGACCTTTACTCCTGCACTGTGACACACGAGATTGACCGC
TACACGGCAATTGCCTATTGGGTACCCAGAACGCCCTGCCTTCGGATCTCCTAGAGAATGCCCTGTGTG
GTGTGGCCTTTGGCCTCGGTGTGCTAGGCACCATCATGGGCATTGCTTCTTCTCTGCTCCAGAGACC
GTGCTCAGGTGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203397 representing NM_010386
 Red=Cloning site Green=Tags(s)

MEHEQKSGAVLLRLLRLLWLLPHSWAVLEASTPVFWDQPNTFRHTLFCQDGPINIGLSETYDEDELFS
 FDFSQNTRVPRLPDFAEWAQGGDASAIAFDKSFCEMLMREVS PKLEGQIPVSRGLPVAEVFTLKP LEFG
 KPNTLVCFISNLFPPTLTVNWQLHSAPVEGASPTSISAVDGLTFQAFSYLNFTPEPFDLYSCTVTHEIDR
 YTAIAYWVPQNALPSDLLLENALCGVAFGLGVLGTIMGIVFFLCSQRPCSGD

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_010386

ORF Size: 783 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_010386.4 , NP_034516.3
RefSeq Size:	1658 bp
RefSeq ORF:	786 bp
Locus ID:	14998
Cytogenetics:	17 17.98 cM