

Product datasheet for **MG203434**

Aqp8 (NM_007474) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Aqp8 (NM_007474) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Aqp8
Synonyms:	A1255744; AQP-8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203434 representing NM_007474 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTGGGGAGACCAATGTGTAGTATGGACCTACCTGAGGTCAAGGTGAAGACCAGCATGGCTGGCA
GATGCCGTGTGTTCTGGTATGAGCAATACGTACAACCATGTATAGTGAACTTGTGGCTCCGCTCTCTT
CATCTTCATCGGATGTCTATCGGTCATTGAGAATAGTCCGAATACTGGGCTCCTGCAGCCTGCCTTGGCT
CATGGCTGGCCTTGGGCTCATCATTGCTACCTTGGGAACATCAGCGTGGACACTTCAACCCTGCTG
TGTCGCTGGCAGTCACAGTGATCGGAGGCTCAAGACCATGCTGCTAATCCCTATTGGATCTCCAGCT
ATTTGGAGGGCTGATTGGGCTGCCTTGGCTAAAGTGGTGGTCCAGAGGAAAGATTCTGGAATGCATCT
GGGGCAGCCTTTGCCATCGTCCAGGAGCAGGAGCAGGTGGCAGAAGCCCTGGGGATAGAGATCATTCTGA
CAATGCTGTTGGTATTGGCTGTGTGTATGGGTGCTGTCAATGAGAAGACAATGGGCCCTTTAGCCCCATT
CTCCATTGGCTTCTGTGCTATTGTGGACATCCTGGCAGGTGGTAGCATCTCTGGAGCCTGCATGAACCT
GCTCGTGCCTTTGGACCTGCTGTGATGGCTGGCTACTGGGACTTCCATTGGATCTACTGGCTGGCCCCC
TCCTGGCTGGCCTCTTTGTAGGACTGCTCATTAGGCTCCTCATCGGAGATGAGAAAACCCGCTAATTCT
AAAGTCGAGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203434 representing NM_007474
 Red=Cloning site Green=Tags(s)

MSGETPMCSDLPEVKVKTSMAGRCRVFWYEQYVQPCIVELVGSALFIFIGLSVIENSPNTGLLQPALA
 HGLALGLIATLGNISGGHFNPAVSLAVTVIGGLKTMLLIPYWISQLFGGLIGAALAKVVSPEERFWNAS
 GAAFAIVQEQQVAEALGIEIILTMLLVLAVCMGAVNEKTMGPLAPFSIGFSVIVDILAGGSISGACMNP
 ARAFGPAVMAGYWDFHWIYWLGPLLAGLFFVGLLIRLLIGDEKTRLILKSR

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_007474

ORF Size: 1483 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_007474.1
RefSeq Size:	1448 bp
RefSeq ORF:	786 bp
Locus ID:	11833
UniProt ID:	P56404
Cytogenetics:	7 F3
Gene Summary:	Forms a water-specific channel; mercury-sensitive. Also permeable to urea but not to glycerol.[UniProtKB/Swiss-Prot Function]